# **SafeMask**<sup>®</sup> TailorMade<sup>™</sup>





| DETAILED MASK DIM | ENSIONS |
|-------------------|---------|
| Height of Earloop | 62,5 mm |
| Pleat Depth       | 15 mm   |



# S Medicom TECHNICAL DATA SHEET

# **FEATURES & BENEFITS**


- Mask and earloops are latex free, safe for wearers with latex allergy or sensitivity.
- Soft, stretchy, ultrasonically sealed earloops prevents irritation and allows easy removal with reduced risk of self-contamination.
- Soft, non-woven material allows for easy breathing.
- Adjustable nose and chin pieces form a seal around the face to reduce gaping.

#### RECOMMENDED USE

• Procedure masks are intended to be worn by personnel during medical procedures to protect both the patient and the operating personnel from transfer of microorganisms, body fluids and particulate material transfer. Latex free.

# PRODUCT DETAILS

Product Name: Medicom SafeMask<sup>®</sup> TailorMade™ **Earloop Face Mask** Product Type: Single-use, non-sterile, three-layer construction Earloop: Polyester and spandex Nosepiece: Aluminum Chin-Wire: Aluminum Colour: Blue Quantity per Inner Carton: 50 China Country of Origin:



| ASTM F2100 - LEVEL 3 BARRIER PROTECTION               |              |   |  |  |  |  |  |
|---|--------------|---|--|--|--|--|--|
| PERFORMANCE FACTORS                                   | TEST METHODS | RESULT                                    |  |  |  |  |  |
| Bacterial filtration efficiency (BFE)                 | ASTM F2101   | ≥ 98%                                     |  |  |  |  |  |
| Sub-micron particulate efficiency at<br>0.1µm (PFE) / | ASTM F2299   | ≥ 98%                                     |  |  |  |  |  |
| Resistant to penetration by synthetic blood           | ASTM F1862   | @ 160 mm Hg                               |  |  |  |  |  |
| Differential Pressure (∆P)                            | EN 14683     | < 6.0 mm H <sub>2</sub> O/cm <sup>2</sup> |  |  |  |  |  |
| Flame Spread  | 16 CF 1610   | Class 1                                   |  |  |  |  |  |

#### NON-WOVEN SPUNBOUND OUTER & INNER LAYER

Provides maximum protection against fluid penetration. Outer and inner layers made from spunbond polypropylene (SBPP).

#### **HIGH-QUALITY MEDIA**

Ensures optimal filtration and breathability. Made from melt blown polypropylene (MBPP).

# S Medicom<sup>®</sup>

#### SOFT, STRONG, ULTRASONICALLY SEALED EARLOOPS

Are attached to the outside of the mask to prevent pulling and irritation and allow easy removal with less risk of self-contamination.

#### **SOFT-BOND ULTRASONICS**

Ensure solid construction and eliminate holes and defects.

#### SHINGLE PLEAT DESIGN

Prevents fluid pooling and cross-contamination. Mask must be worn with pleats facing down.

#### ADJUSTABLE, ENCLOSED NOSEPIECE

Forms a strong seal to provide customized protection and comfort.

# CERTIFICATION & STANDARDS

ISO 13485 certified manufacturer. Class 1 Medical Device. Health Canada registered medical device manufacturer.

## STORAGE CONDITIONS

These masks should be disposed of after each use. Store in a cool, dry place. Avoid excessive heat (40 °C or 104 °F).



## DONNING INSTRUCTIONS

- Fan to open pleats.
- Shape before applying to ensure a better fit.
- o Pull earloops over each ear.
- With two fingers on either side of your nose, firmly shape around your nose and across your cheeks, concentrating on the bridge of your nose, to secure the seal.
- Pinch the nosepiece and make sure the base of the mask sits comfortably below your chin.

## DOFFING INSTRUCTIONS

- o Pull the earloops off your ears.
- Pull the mask away from your face, taking care not to touch its outer surface.
- Still holding only the earloops and taking care not to touch outer surface, dispose of mask in a clinical waste bin.
- o Always wash your hands thoroughly after removing a mask.



| Appendix 1 (List of items) |  |   |       |   |             |                           |                           |  |  |  |
|----------------------------|--|---|-------|---|-------------|---------------------------|---------------------------|--|--|--|
| Item Code                  | Description  | Packaging<br>Information                          | Color | Specific<br>Attribute   | Treatment   | Marketed in<br>USA as     | Marketed in<br>Canada as  |  |  |  |
| 2072                       | Medicom SafeMask Tailor<br>Made Earloop Face Masks<br>with Chin Wire | 50 Masks/Box, 10<br>Boxes/Case, 500<br>Masks/Case | Blue  | Aluminum Nose-piece,<br>Flat Knitted Earloop,<br>Regular Earloop Length,<br>Chin Wire | Non-Sterile | Medical Device<br>Class 2 | Medical<br>Device Class 1 |  |  |  |

AMD Medicom Inc. | 2555 Chemin de l'Aviation, Montreal, Quebec, Canada H9P 2Z2