


# INSTRUCTIONS FOR USE

Personal Protection



## Surgical Earloop Face Masks

How to Properly Wear Your Crosstex Face Mask:



Step 4 for SecureFit™ Masks only

Ref.#'s GCLBL, GCLPK, GCBL\*, GCPK\*, GCAPTK\*, GCIBLSF, GCILVSF, GCIPKSF, GCIBL, GCIGR, GCILV, GCIPK, GCISA, GCITE, GCITZ, GCITQ, GCIBL100, GCIBL300, GCICXB, GCICXS, GCICXZ, GCICXT, GCIPWB, GCIPWS, GCIPWZ, GCIPWT, GCICXBSF, GCIPWBSF, GPLUSBL, GPLUSKA, GPLUSPK, GPLUSWHSF, GCPBLSF, GCPLVSF, GCPPKSF, GCPBL, GCPLV, GCPPK, GCPYE, GCPBL100, GCABL\*, GCAPK\*, GCFCXSSF, GCFCXSFSF, GCPWSSF, GCFCXUSF, GCFCXS, GCFCXSF, GCPWS, GCFCXU, GCFCX, GCPW

### Intended Use:

Crosstex™ Surgical Masks are intended for use in infection control practices to minimize contamination caused by inhaled and exhaled microorganisms and reduce the potential exposure of the wearer to blood and body fluids. This is a single use, disposable device and provided non-sterile.

### Indications for Use:

Crosstex™ Surgical Masks are intended for use in infection control practices to minimize contamination caused by inhaled and exhaled microorganisms and reduce the potential exposure of the wearer to blood and body fluids. This is a single use, disposable device and provided non-sterile.

### Contra-Indications:

Do not re-use, this is a single-use disposable device.  
Do not use near flammable anesthesia gases.

### Compliance:

US-FDA 510K cleared class II medical device Surgical Mask.  
US-FDA Regulation: 21 cfr 878.4040.  
Meets fluid barrier protection per ASTM F2100 standards.  
Meets applicable ISO 10993 Biocompatibility testing requirements.  
Meets Class 1 Flammability Spread per 16 cfr 1610.  
Devices are manufactured in compliance with 21 cfr Part 820 and ISO 1385:2016 quality systems.



Crosstex International, Inc.  
6789 W. Henrietta Road, Rush, NY 14543 USA  
585.359.0130 | 631.582.6777 | crosstex.com



MDSS GmbH  
Schiffgraben 41  
30175 Hannover, Germany



Do not reuse.














Not made with natural rubber latex.

Made in USA unless otherwise noted. \*Made outside of the USA.  
EARLOOPMASK-IFU / REV D - 07/2021

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.  
©2021 Crosstex International, Inc. All rights reserved.

# CROSSTEX

## Applicable & Definitions Symbols Use:

 crosstex.com	Consult electronic instructions for use
	Catalog number
	Single use only Do not re-use
	Non-sterile
	Manufacturer
	Unique Device Identifier
	Manufacturing lot number
	Date of manufacture
	Medical Device
	Approved for European Union, notified body
	European Union Representative