
## Rectangular P.T. Ball Rebounder 10-3110



- 1. Slide the rebounder handle post through the trampoline frame bracket (insert the post from the front side of the trampoline frame).
- 2. Align the handle post with the post bracket and bolt the two together.
- 3. Mount the ball rack onto the base. The rack can be mounted as shown or on the other side of the post.
- 4.Before using, position the frame bracket at the desired height and angle and tighten the adjustment knob at that height.

## GANDO X Ball Rebounder

- · Great for upper body fitness program with P.T. Balls
- Adjustable angle (0, 28, 33, 39 and 42 degrees) provides variation in training
- Returns the ball at angle, speed and force consistent with throw
- Use to develop agility and range-of-motion



**REF 10-3110** Square rebounder

**REF 10-3132** Rebounder with 5 PT balls (2, 4, 7, 11,15 lb)



REF 10-3112 Round rebounder

**REF 10-3136** Rebounder with 5 PT balls (2, 4, 7, 11,15 lb)

## Suggested CanDo<sup>®</sup> P.T. Balls for Ball Rebounder

<b>SKU</b>	<b>weight</b>	diameter	<b>color</b>
10-3170	1 lb	5"	tan
10-3171	2 lb	5"	yellow
10-3172	4 lb	5"	red
10-3173	7 lb	7"	green
10-3174	11 lb	7"	blue
10-3175	15 lb	9"	black
10-3177	20 lb	9"	silver
10-3178	30 lb	9"	gold



Exercise routines should be determined by a doctor, therapist, coach, or trainer.



Fabrication Enterprises Inc. 250 Clearbrook Rd, Suite 240 Elmsford, NY 10523 (USA) tel: +1-914-345-9300 • 800-431-2830 fax: +1-914-345-9800 • 800-634-5370 FabEnt.com EC REP

AJW Technology Consulting GmbH Königsallee 106 40215 Düsseldorf (Germany)