
Manufactured for:
Dynarex Corporation
10 Glenshaw Street
Orangeburg, NY 10962
USA • www.dynarex.com
Made in China



## Instructions for use

The DynaLift™ Transport Unit is for use by professional staff for lifting and transporting/transferring patients. The DynaLift™ Transport Unit is rated to transport up to 1500 pounds (680 kg). The unit must be placed so handle support straps are to the side opposite the patient. It is recommended to use at least four people to provide safe transport/transfer, thus reducing the possibility of accident or injury. Grab handles in line with the people lifting on opposite sides of the unit – **DO NOT** "Cross-Handle". The DynaLift™ Transport Unit is not designed for use without a backboard when there is a need for spinal stabilization or traction of the patient being transported. Reinforcement straps must face downward!



The DynaLift™ Transport Unit is intended for limited use. If the unit is damaged, frayed, cut, or soiled, replace with a new unit.

- Do Not use if the unit is punctured, torn, frayed, or excessively worn.
- Do Not machine wash or dry.
- Avoid contact with sharp objects.
- Clean soiled unit with a damp cloth. Soap, detergent, or disinfectant can be used.
- · Only use handles for lifting.
- Use care to ensure grip is secure prior to lifting.
- Do Not store in contact with heat sources greater than 200°F (93.3°C).
- Avoid dragging over rough surfaces.
- Use a minimum of four people to lift patient.
- . Do Not attach to mechanical lifts.