

STEP 5

Let your Power Seat lower you into the armchair or sofa.

Please be sure to sit in the center of your Power Seat cushion so that you will be lowered smoothly.



Now sit down on your Power Seat keeping both feet on the floor and holding on to the arm(s).



PUSH DOWN gently on the Power Lever to lower.

Power Seat will stop automatically when it is fully lowered.

Tip: If you are finding it hard to sit down straight keeping both feet on the floor, the Power Seat may be lifted too high. Simply lower it a bit and then be seated.

Please read the *Users Guide* for more detailed information on the use and care of your Power Seat.


SAFETY REMINDERS



Never jump up on your Power Seat. Lower the cushion before use if it is raised too high for you.



Do not sit crooked on your Power Seat. Always sit in the center of the cushion, with both feet on the floor.



Keep your Power Seat out of the reach of children.



Do not use near water or outdoors where the cord could get wet.

Power Seat is not recommended for use in motor vehicles, rockers, wheelchairs or other types of chairs with wheels.

READ THIS FIRST




Thank you for choosing Power Seat. This is a mechanical device, so please take a short time to get used to your Power Seat and practice using it.

Also, please read the *Users Guide* thoroughly for more detail and easy-care instructions.

STEP 1

Choose an armchair or the end of a sofa for your Power Seat.



You should only use your Power Seat in an armchair or at the end of a sofa next to one of the arms as shown above.

STEP 2

Snap on the Power Lever.

Please do not plug your Power Seat in until after you have installed the Power Lever.



You only need one lever to operate the Power Seat. Choose the side that suits you best, or that works best in the armchair or sofa you are using.



Snap the Power Lever onto the metal prong located at either front corner of the cushion as shown.



Angle the levers towards the front of the Power Seat, as shown in the diagram.

This angle makes it easier to find the correct fit for the lever.

Listen for a click so you know the lever has snapped into place.

STEP 3

Learn the 3 lever positions: UP, DOWN and STOP

Please do not sit on the Power Seat until after you have learned the three lever positions.

You may now plug in your Power Seat




PULL UP gently on the lever to raise the Power Seat.

Release the lever. No need to hold on to it.



PUSH DOWN on the lever to lower the Power Seat.

Release the lever.



To STOP, move the lever to the middle position.

Practice until you are very comfortable using all three positions of the Power Lever before going on to Step 4.

STEP 4

Sit on your Power Seat and let it lift you out of the armchair or sofa.

Start with your Power Seat lowered fully before being seated.



Be sure to sit in the center of the cushion and keep your feet slightly forward so that you will be lifted straight up.

PULL UP on the Power Lever to rise.



Hold onto an arm of the chair or sofa during use. Also keep both feet flat on the floor at all times.



Your Power Seat will automatically stop when it has reached its full height, or you can stop it at anytime.

Tip: If you have short legs, simply stop your Power Seat at a height that is comfortable for you