

INSTRUCTIONS—DORSAL RESTING SPLINT

INDICATIONS

C.V.A.—Spinal Cord Industries—C.P.

APPLICATION

Open straps.

With forearm section on top, pass fingers through opening between hand and forearm sections. Fingers on top, thumb on bottom, and web space against hand section.

Position finger separators. Secure hand strap and position strap pad.

Secure forearm strap and position pad. Pad can be removed and trimmed.



Tightening forearm will bring forearm section and wrist into position.

EVALUATE FIT, FUNCTION, AND INSTRUCT PATIENT ON USE.

SIZING CHART

Size	Width at M.P. Joints
Pediatric	2 1/4" - 2 3/4"
Small	2-3/4" - 3-1/4"
Medium	3-1/4" - 3-3/4"
Large	3-3/4" and up

* Replacement Liners Available



ADJUSTMENTS

The lateral portion of the forearm section and hand section are malleable aluminum and can easily be shaped. The wire frame connecting the hand and forearm section is strong and requires more force to bend.

Additional padding—hook and loop compatible (provided)—place between liner and splint.

Straps

The straps are secured to splint with hook material and can be repositioned to accommodate use by right or left hand.

Forearm strap length adjustment—Remove strap from hook “sandwich”. Cut end of strap WITHOUT hook to length and replace in “sandwich”.

Optional Wrist Strap—Attach adhesive backed hook to splint in desired position. Attach strap, adjust length and position pad.

CLEANING

Liner—Remove from splint. Hand wash in cold water with mild detergent. Rinse thoroughly and air dry completely before reusing.

Splint—Wipe with damp towel and mild detergent. Towel dry completely before reusing.

Avoid placing splint in hot environment—next to heater—hot car

Immediately contact your healthcare professional if you are experiencing skin irritations, pain, loss of sensations or if you have any questions or concerns.

Single patient use. To be used only when ordered by a physician and applied by qualified medical professional



Authorized Representative: Medica Surgical Innovations Limited, Moorgate Street, Blackburn BB2 4PB, Lancs UK