DxH 500 Series Calibrator

REF B36880

TABLE OF CONTENTS

English	2
Français (FR)	5
Deutsch (DE)	8
Italiano (IT)	11
Español (ES)	14
Português (Portugal)	17
Dansk (DA)	20
Svenska (SV)	23
Norsk (NO)	26
Ελληνικά (EL)	29
日本語 (JP)	32
中文 ZH-CH	35
Lietuviškai (LT)	38
Magyar (HU)	41
Polski (PL)	44
Čeština (CZ)	47
Slovenčina (SK)	50
한국어 (KO)	53
Türkçe (TR)	56
Русский (RU)	59
Hrvatski (HR)	62
Български (BG)	65
中文 ZH-TW	68
Română (RO)	71
Србија (SR)	74
Українська	77
Português (Brasil)	80
Nederlands	83
Tiếng Việt (VI)	86

DxH 500 Series Calibrator

REF B36880

**For *In Vitro* Diagnostic Use
Rx Only in the U.S.A.**

INTENDED USE

DxH 500 Series Calibrator is designed for the determination of calibration factors for the DxH 500 Series Systems in conjunction with specific DxH 500 Series reagents. Refer to the Table of Expected Results for specific instrument models.

SUMMARY AND PRINCIPLE

Hematology analyzers require periodic calibration in order to generate accurate patient results. This calibrator is a stable, whole blood preparation that can be used to verify and adjust calibration of select hematology instruments.

Calibrator values for DxH 500 Series Calibrator are derived from replicate testing on instruments operated and maintained according to the manufacturer's instructions. Instruments utilized for assay value assignment are whole blood calibrated using values determined by reference methods, as detailed in the Reference Methods Section of this IFU.

REAGENTS

DxH 500 Series Calibrator consists of human erythrocytes, mammalian leukocytes and mammalian platelets suspended in a plasma-like fluid with preservatives.

WARNING AND PRECAUTIONS

POTENTIAL BIOHAZARDOUS MATERIAL. For *in vitro* diagnostic use. Each human donor/unit used in the preparation of this product has been tested by a FDA licensed method/test and found to be negative or non-reactive for the presence of HBsAg, Anti-HCV, NAT testing for HIV-1, HCV (RNA) and HIV-1/2. Each unit is also negative by a serological test for Syphilis (RPR or STS). Because no test method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious. When handling or disposing of vials follow precautions for patient specimens as specified in the OSHA Bloodborne Pathogen Rule (29 CFR Part 1910, 1030) or other equivalent biosafety procedures. DxH 500 Series Calibrator is intended for use only by trained personnel.

GHS HAZARD CLASSIFICATION		
DxH 500 Series Calibrator	WARNING 	Causes skin irritation.
		Causes serious eye irritation.
		Wear protective gloves, protective clothing and eye/face protection.
		IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
		If skin irritation occurs: Get medical advice/attention.
		If eye irritation persists: Get medical advice/attention.
Potassium Hydroxide <1%		

STORAGE, STABILITY, AND DISPOSAL

Store DxH 500 Series Calibrator upright at 2-8°C (35-46°F) when not in use. **Protect tubes from overheating and freezing.** Unopened tubes are stable until the expiration date. Opened tubes are stable for 5 days, provided they are handled properly. Dispose of waste product, unused product and contaminated packaging in compliance with federal, state and local regulations.

EVIDENCE OF DETERIORATION

After mixing, product should be similar in appearance to fresh whole blood. In unmixed tubes, the supernatant may appear cloudy and reddish; this is normal and does not indicate deterioration. Other discoloration, very dark red supernatant or unacceptable results may indicate deterioration. **Do not use the product if deterioration is suspected.**

INSTRUCTIONS FOR USE

A. Mixing and handling directions:

1. Remove tube from the refrigerator and verify that the lot number on the tube matches the lot number on the Table of Expected Results.
2. Allow the tube to warm at room temperature (15-30°C or 59-86°F) for 15 minutes before mixing.
3. To mix, hold a tube horizontally between the palms of the hands. **Do not pre-mix on a mechanical mixer.**
 - A. Roll the tube back and forth for 20-30 seconds; occasionally invert the tube.
 - B. Continue to mix in this manner until the red cells are completely suspended. Tubes stored for a long time may require extra mixing.
 - C. Gently invert the tube 8-10 times immediately before sampling.

B. Analyze Calibrator:

1. Analyze calibrator according to the procedure in the Instructions for Use for your instrument.
2. Prime the instrument once by aspirating the calibrator sample. The first run should be excluded as a prime.
3. Clean residual material from the cap and rim if the tube is opened for sampling. Replace the cap tightly.
4. Compare the mean value for each parameter to the assigned value.
 - A. If the difference is within the Acceptable Range, calibration is optional.
 - B. If the difference exceeds the Acceptable Range, calibration may be needed.
5. Acceptable ranges given on the assay sheet are intended as guidelines, but not absolute limits, for evaluating instrument calibration. Acceptable calibration should be established by each laboratory.

C. Adjust instrument calibration and verify results:

1. Calibrate the instrument by using the calibration adjustment procedures described in the Instructions for Use for your instrument.
2. Verify calibration by analyzing a new tube of calibrator and repeat section A & B.

LIMITATIONS

The performance of this product is assured only if it is properly stored and used as described in this insert. Incomplete mixing of a tube prior to use invalidates both the sample withdrawn and any remaining material in the tube.

PRODUCT AVAILABILITY

DxH 500 Series Calibrator

REF B36880 - 2 x 2 mL

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

For additional information or if damaged product is received, call Beckman Coulter Customer Service at 800-526-7694 (USA or Canada) or contact your local Beckman Coulter Representative.

ADDITIONAL INFORMATION

Glossary of Symbols is available at techdocs.beckmancoulter.com (PN C07163)

May be covered by one or more pat. - see www.beckmancoulter.com/patents

SDS Safety Data Sheet is available at techdocs.beckmancoulter.com

Revision History

Revision AF, 01/2018

Changes were made to:

- Added new Serbian language
- Added new Vietnamese language

Revision AG, 09/2018

Changes were made to:

- Update Russian language

Revision AH, 09/2018

Changes were made to:

- Added new Patent Information Statement
- Added new Dutch language

REFERENCE METHODS

1. **WBC:** A series of 1:500 dilutions are made with calibrated glassware. Counting is performed on a Coulter Counter Z series instrument. All counts are corrected for coincidence.
 2. **RBC:** A series of 1:50,000 dilutions are made with calibrated glassware. Counting is performed on a Coulter Counter Z series instrument. All counts are corrected for coincidence.
 3. **HGB:** Hemoglobin value is determined by spectrophotometric procedure according to CLSI Standard H15-A3 and is traceable to ICSH/WHO International Haemiglobincyanide Standard.
 4. **HCT:** Packed cell volume (PCV) is measured by the microhematocrit procedure according to CLSI Standard H7-A3. No correction is made for trapped plasma.
 5. **MCV:** Calculated using the reference RBC and HCT.
 6. **PLT:** A series of 1:126 dilutions are made using calibrated glassware in 1% ammonium oxalate. Platelets are counted using a hemocytometer and phase contrast microscopy.
-

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Bản sửa đổi AG, 09/2018
Beckman Coulter Ireland Inc., Lismeehan, O'Callaghan's Mill's Co. Clare, Ireland (001) 703-527-3887
Manufactured for Beckman Coulter, Inc.