

# ENZYCLEAN®

Protease Enzyme  
Low Suds Detergent



Enzyclean® Protease Enzyme Low Suds Detergent is specially formulated for use as a presoak or manual detergent and for use in ultrasonic washers, automatic washing equipment and evacuators. Enzyclean® is safe for glass, plastic, rubber and all types of metal surgical instruments, as well as flexible and rigid fiber optic equipment.

- **Protease enzymes effectively break down proteins.**
- **Compatible with all equipment & pump systems.**
- **Product safe and compatible with all metals.**

## Protease Enzyme Low Suds Detergent

Order No.	Size	Physical Properties
EBL1	1 Gallon (3.8 Liters)	<b>Form:</b> Liquid <b>Color:</b> Green <b>Odor:</b> Fresh <b>pH Concentration:</b> 8.7

Uses	Temperature
<input checked="" type="checkbox"/> Presoak <input checked="" type="checkbox"/> Manual <input checked="" type="checkbox"/> Ultrasonic <input checked="" type="checkbox"/> Evacuators	<b>Usage:</b> 60°-140° F (16°-60° C) <b>Storage:</b> Cool, dry area
<input checked="" type="checkbox"/> Automatic Washer <input checked="" type="checkbox"/> Automatic Endoscope Reprocessor <input type="checkbox"/> Cart Washer	

Product Features	Consumer Benefits
<b>Protease Enzymes:</b>	Protease enzymes effectively break down proteins.
<b>Rust Inhibitors:</b>	Reduces likelihood of rust, tarnish or corrosion.
<b>Rinsing Agents:</b>	Rinses spot & residue-free.
<b>Highly Concentrated:</b>	Low cost per use.
<b>Neutral pH:</b>	Product safe and compatible with all metals.
<b>Low Sudsing:</b>	Compatible with most equipment & pump systems.
<b>Bacteriostatic:</b>	If stored properly, not favorable to bacterial growth.
<b>Environmentally-Friendly:</b>	Biodegradable, phosphate-free, VOC free formula.

Procedures	Directions for Use
<b>Directions:</b>	Wear appropriate personal protective equipment based on safe work practices.
<b>Presoaking, Manual and Ultrasonic Cleaning:</b>	Add 1/8 fl. oz. (3.7 mL) to 2 fl. oz. (59 mL) per gallon (3.8 L) of warm water. Allow to soak for 2-3 minutes or longer. Discard solution after each use.
<b>Automatic Washing Equipment:</b>	Add 1/8 fl. oz. (3.7 mL) to 2 fl. oz. (59 mL) per gallon (3.8 L) of water. Recommended temperature range: 90°-140° F (32°-60° C).
<b>Evacuators:</b>	Add 1/8 fl. oz. (3.7 mL) to 2 fl. oz. (59 mL) per gallon (3.8 L) of water. Suggested maintenance twice weekly.
	<b>Dilutions vary depending on temperature, water quality and equipment. Validated cleaning is recommended.</b>