

Safe-Flo™ Saliva Ejectors and Valves



Protect your patients from evacuation line backflow

Evacuation line health risk is real.

Backflow in low-volume evacuation lines can occur when the pressure in the patient's mouth is less than in the evacuator. This condition can result under several commonly occurring situations¹:

- Patient creates a seal around the saliva ejector by closing their lips around the tip
- Saliva ejector tip becomes occluded from cheek, tongue or other blockage
- Suction tubing is positioned above the patient's mouth
- Simultaneous use of other evacuation (HV) equipment

Patented Backflow Prevention²


The unique, one-way valve provides a physical barrier to prevent backflow from occurring.

Patient Safety


Keep all your patients safe from the potential risk of contaminated backflow. Provides added protection for geriatric patients who may experience a variety of oral soft tissue changes, as well as pediatric and special needs patients who may inadvertently close their lips around the saliva ejector.

Seamless Workflow

The integrity of suction and performance of the saliva ejector or HVE tip is not affected, so there is no need to change how you work in order to protect from backflow.



Unique, **one-way valve** immediately closes if vacuum pressure changes, creating a physical barrier to prevent backflow.




Available as a complete saliva ejector/valve system or sold separately for use with your preferred brand of saliva ejectors or HVE tips.

DID YOU KNOW?

Studies show contaminated **backflow could potentially occur 21% of the time³**. Safe-Flo™ Saliva Ejectors and Valves provide backflow prevention you can rely on to protect all of your patients.

Backflow may contain:



BACTERIA



BLOOD



SALIVA



TISSUE

	Description	Quantity	Part Code
A	Safe-Flo™ Saliva Ejector Clear/Blue	100/Bag	M1001CLBU
B	Safe-Flo™ Saliva Ejector, White/Blue	100/Bag	M1000WHBU
C	Safe-Flo™ SE Valve, Blue	100/Bag	M1010SUSE
D	Safe-Flo™ HV Valve, Gray	100/Bag	M1030HVSE
E	Safe-Flo™ Adapter, Blue/Gray	50/Bag	M1020ADAP

Crosstex is now a proud member of

HuFriedyGroup
The Best In Practice

Safe-Flo™ Saliva Ejectors and Valves




Safe-Flo™ Saliva Ejectors

- For backflow prevention in low-volume evacuation lines.
- Fully assembled, all-in-one saliva ejector with Safe-Flo™ Valve replaces your current saliva ejector brand.
- Slightly shortens the exposed length of a saliva ejector, making it less flexible, allowing for easier tissue retraction and better access.




Safe-Flo™ Saliva Ejector Valve

- For backflow prevention in low-volume evacuation lines.
- Use with your preferred saliva ejector brand.




Safe-Flo™ High Volume Evacuation Valve

- For backflow prevention in high-volume evacuation (HVE) lines.
- Use with your preferred HVE tip brand.



Safe-Flo™ Adapter

- Allows you to use the HVE line with your preferred saliva ejector brand for backflow prevention.
- Provides a higher level of aspiration to the saliva ejector.



Use a Safe-Flo™ Valve with any evacuation line.

Low Volume



Saliva Ejector Tip



Safe-Flo™ SE Valve



On/Off Switch



Evacuation Line Connector

High Volume



HVE Tip



Safe-Flo™ HV Valve



On/Off Switch



Evacuation Line Connector

Safe-Flo™ Valves are single-use disposable, one-way valves that stop backflow from evacuation lines, preventing cross-contamination.

Call your preferred distributor to order, or learn more at [HuFriedyGroup.com/Safe-Flo](https://www.HuFriedyGroup.com/Safe-Flo)

Not made with natural rubber latex.

1) Center for Disease Control and Prevention. Guidelines for Infection Control in Dental Health Care Settings-2003. MMWR 2003;52 (no. RR-17); [31].

2) US Pat. No. 8,256,464, 8,714,200, 8,752,579 and foreign counterparts.

3) Miller C. Back flow in low-volume suction lines may lead to potential cross-contamination. RDH. 1996; 16(1):30.

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | [HuFriedyGroup.com](https://www.HuFriedyGroup.com)

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, their affiliates or related companies, unless otherwise noted.

Marks not registered in all jurisdictions.

©2023 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-807 / 0523



Crosstex is now a proud member of
HuFriedyGroup
The Best In Practice