



 THERABAND®

 THERABAND®  
FLEXBAR® | PerformanceHealth.com

# FlexBar® EXERCISE AND TREATMENT GUIDE



Proven Science, Trusted Performance.



## What is FlexBar®?

FlexBar® is a lightweight, easy-to-grip, portable resistance exerciser. It is engineered for wrist, forearm and shoulder strengthening and recovery.

DESIGNED & MADE IN THE  
**USA**




## How can I use FlexBar® in my practice?

- **Strengthening:** improves grip, wrist, shoulder and thumb strength
- **Oscillation:** allows oscillation movements for the elbow and shoulder for neuromuscular training and balance training
- **Mobilization:** provides soft tissue mobilization and joint mobilization for the hand, neck, foot, fingers and toes



## What type of material is FlexBar®?

FlexBar® is made from dry natural rubber. Each bar is 12 inches long, with resistance levels that increase with each of **four color-coded diameters.**



Represents typical values.

# FlexBar® – EXERCISE AND TREATMENT

**Proven Performance:** Clinicians worldwide utilize FlexBar® for the following rehabilitation exercises.

## TYLER TWIST

### Tennis elbow

**Treatment:** *Lateral epicondylitis*

Eccentric contraction of wrist extensors to increase forearm and wrist strength for lateral elbow pain treatment.



STEP

1

Hold FlexBar® upright with hand of injured arm, wrist extended.



STEP

2

While holding bar, grasp FlexBar® as shown with opposite hand.



STEP

3

Twist FlexBar® by flexing non-injured wrist.



STEP

4

Bring arms forward with elbows extended while maintaining twist in FlexBar®.



UNTWIST



STEP

5

Slowly untwist FlexBar® by allowing injured wrist to move as it unwinds.

# REVERSE TYLER TWIST

## Golfers elbow

**Treatment:** *Medial epicondylitis*

Eccentric contraction of the wrist flexors and pronators for medial elbow pain treatment.



STEP  
1

Bend injured elbow and hold FlexBar® parallel to the ground.



STEP  
2

Lift opposite elbow upward and rotate forearm so palm faces away from body.



STEP  
3

Grasp FlexBar® with non-injured hand facing away from the body and downward.



STEP  
4

Twist FlexBar® with non-injured arm as you stabilize the injured arm.



STEP  
5

Hold both wrists steady as you extend both elbows in front of you. The wrist on your injured side should be flexed toward you, and the other wrist extended.



STEP  
6

Slowly untwist the FlexBar® with the injured side while maintaining tension with the non-injured side.

# SHOULDER OSCILLATION

## Shoulder rehabilitation

**Treatment:** Increase muscle activity and restore muscle balance in different shoulder positions and planes of oscillation.



STEP

1

Stand with feet shoulder-width apart, arms alongside body. Back and trunk remain straight with belly button in tight. Knees are "loose" – not hyperextended.



STEP

2

With one hand, grasp the end of a FlexBar® and lift that arm into flexion in front of body, with hand grip at the base of the FlexBar® and arm just below shoulder height.



STEP

3

Oscillate FlexBar® for 20-30 seconds. Rest and repeat.

# WRIST FLEXION INLINE

## Wrist pain

**Treatment:** Strengthen the wrist and finger flexors (forearm).



STEP

1

Grasp FlexBar® vertically with both hands, elbows bent in front of body, palms facing trunk. Both hands squeeze FlexBar® to hold a tight grip.



STEP

2

With hand of wrist to be exercised (top one in picture), push against FlexBar® to curl palm in toward elbow. Other fist stabilizes against the motion. Hold for 2-3 seconds. Slowly return to starting position (neutral). Repeat.

# Evidence-based results

**THERABAND® products are backed by more than 800 evidence-based research articles** and references with quantifiable, proven results.

Below are examples of FlexBar®-specific research.



**Muscle activity of upper extremity increases during oscillation exercise using FlexBar®**

Page P et al., 2004

**Summary:** Using surface electromyography, this study investigated the muscle activation of the upper extremity during a THERABAND® FlexBar® oscillatory exercise in different shoulder positions and planes of oscillation.

**Results:** The scapular muscles exhibited the highest activation in the scaption/sagittal condition, suggesting that this condition be used in the rehabilitation of scapular stabilizers. The ratio between lower trap and upper trap activation demonstrated values of 1.6 to 2.3, indicating the flexion position may be optimal for improving this ratio.



**Breakthrough strength training with FlexBar® may help treat tennis elbow**

Tyler TF et al., 2010

**Summary:** In the study, one group of patients with tennis elbow received traditional physical therapy, while the other received traditional physical therapy with the addition of the “Tyler Twist” eccentric exercise using THERABAND® FlexBar®.

**Results:** The patient group utilizing the Tyler Twist exercise in their treatment plan reported an 81% improvement in elbow pain and a 72% improvement in strength, while the control group showed little or no improvement.



**Treating medial elbow pain diagnosed as golfer's elbow using FlexBar®**

Hogan, D et al., 2014

**Summary:** Twenty patients with medial elbow pain diagnosed as golfers elbow were followed in this study. All the patients had failed previous conservative treatments such as medication, injections and physical therapy.

**Results:** A FlexBar® exercise prescription of three sets of 15 repetitions daily for approximately six weeks appeared to be an effective treatment in the majority of patients who had already failed a previous intervention for this disorder.

“The use of the THERABAND® FlexBar® in concert with a physical therapy treatment regimen for wrist extensor strengthening is effective in reducing pain and improving wrist and forearm strength. We use the FlexBar® regularly in our clinic and for the treatment and training of elite level tennis players at tennis facilities around the world.”

**- Todd S. Ellenbecker, PT, DPT, MS, CSCS, FAPTA**  
Director of Shoulder Rehabilitation & Director of Clinical Research and Publication



**Scan the QR code to view these and other articles.**

Caution: This product contains natural rubber latex, which may cause allergic reactions. Always read the label. Follow the instructions for use. THERABAND®, FlexBar® and the Color Pyramid & Associated Colors are all trademarks of Performance Health and/or its subsidiaries and may be registered in the United States and other countries. ©2023 All rights reserved. Unauthorized use is strictly prohibited. P11819-R00

 **THERABAND®**

THERABAND.COM