6000: FRACTURE KIT: INSTRUCTIONS FOR USE

SAVE THESE INSTRUCTIONS FOR FUTURE USE

Note: The most current version of these instructions can be found online at www.grahamfield.com.

- Λ WARNING: Important! Read and understand these instructions before using the 6000 Fracture Kit. This Fracture Kit should be used only by a person trained in providing emergency medical attention. If the 6000 Fracture Kit is not properly applied, personal injury and/or damage to the Fracture Kit could result.
- Use only Grafco replacement parts. Non-Grafco replacement parts could cause personal injury and/or damage to the Fracture Kit.
- assembly or use of this product.

FRACTURE KIT COMPONENTS (SEE PICTURE AT RIGHT)

Instructions for Use, Splints:

- 1. Select appropriate splint.
- 2. Support the injured limb above and below the suspected break.
- 3. Position the splint along the long axis of the limb.
- 4. Ensure there is a pulse before attaching hook & loop straps.
- 5. Attach hook & loop straps and support the limb for transport.
- 6. Recheck to ensure pulse on injured limb.
- 7. Transport patient to medical care.

Instructions for Use, Cervical Collar:

- 1. While supporting patient's head, place the padded hook & loop strip behind patient's neck.
- 2. Gently slide the curved chin piece firmly in place at front of patient's neck.
- Attach the loop strap to the chin piece hook strap. 3.
- 4. While maintaining head support, secure the patient to the transport device. Use sandbags or other immobilizers as needed. Arrange for transport to medical care.

LIMITED WARRANTY

SCOPE OF WARRANTY GF Health Products, Inc. ("GF") warrants to the original purchaser only that it will replace or repair components, at GF's sole discretion, that are defective in material or workmanship under normal use and GF Health Products, Inc. ("GF) warrants to the original purchaser only that it will replace or repair components, at GF sole discretion, that are derective in material or workmarship under normal use and maintenance. The property warrants are conditioned upon the proper use of the products strictly in accordance with good commercial practice and applicable GF instructions and manuals, including proper use and maintenance. To the extent that a third party warrants a component, GF conveys all of its rights under that warranty to the original purchaser, to the extent permitted. This limited warranty shall only apply to defects that are reported in accordance with the provisions set forth in this warranty document, within the applicable warranty period and which, upon examination by GF or its authorized representative, prove to be a warranty item (See Obtaining Warranty Service below) This limited warranty is not transferable. Labor is not included in the warranty period is as designated above. If a part is replaced under warranty, the original warranty period will not be affected. All other replacement parts will be subject to the warranty period specified. The applicable warranty period shall commence from date of shipment to the original customer, unless there is an expiration date on the component in which case the warranty shall expire on the earlier of warranty period or the expiration date.

OBTAINING WARRANTY SERVICE Customers located in the United States who wish to report a warranty issue, must contact GF directly by calling 1770.368.4700 or by e-mailing a request to cs@grahamfield.com. Customers located outside the United States must contact the Distributor from whom they purchased the products. In both cases, further directions will be provided once the initial contact is made. This limited warranty shall only apply to defects that are reported within the applicable warranty period. Failure to abide by the specific directions will result in denial of the warranty claim. The warranty does not cover and GF shall not be liable for the following:

- Defects, damage, or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, tampering or failure to seek and obtain repair or replacement in a timely manner; 1)
- 2) Products which are not installed, used, or properly cleaned and maintained as required in the official manual for the applicable product;
- Products considered to be of a non-durable nature including, but not limited to: filters, fuses, gaskets, lubricants, and charts; 3)
- 4) Accessories or parts not provided by GF;
- Matching of color, grain or texture except to commercially acceptable standards; Changes in color caused by natural or artificial light: 5) 6)
- Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing, in advance, by GF; 7)
- 8) Any labor or shipping charges incurred in the replacement part installation or repair;
- Costs and expenses of regular maintenance and cleaning; and 10) Representations and warranties made by any person or entity other than GF.

ENTIRE WARRANTY, EXCLUSIVE REMEDY AND CONSEQUENTIAL DAMAGES DISCLAIMER

ENTIFE WARRANTY IS GFS ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. GF MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTY IS GFS ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. GF MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTY IS OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IF ANY MODEL OR SAMPLE WAS SHOWN TO THE CUSTOMER, SUCH MODEL OR SAMPLE WAS USED MERELY TO ILLUSTRATE THE GENERAL TYPE AND QUALITY OF THE PRODUCT AND NOT TO REPRESENT THAT THE PRODUCT WOULD NECESSARILY CONFORM TO THE MODEL OR SAMPLE IN ALL RESPECTS. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF THE DEFECTIVE PARTS. GF SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO: DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER, OR EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. CERTAIN STATES MAY CONFER ADDITIONAL RIGHTS REGARDING WARRANTIES AND IN THOSE STATES GFS LIABILITY AND THE LIABILITY OF GFS SUDDIERS. SHALL BE LIMITED TO THE FEILUREST PROTOCHTED RY LAW SUPPLIERS, SHALL BE LIMITED TO THE FULLEST EXTENT PERMITTED BY LAW.

The warranties contained herein, together with GF's current Terms and Conditions, contain all the representations and warranties with respect to the subject matter of this document, and supersede all prior negotiations, agreements and understandings with respect thereto. The recipient of this document hereby acknowledges and represents that it has not relied on any representation, assertion, guarantee, warranty, collateral contract or other assurance, except those set out in this document.

For additional information on this product or this warranty, please contact a GF Customer Service Representative

NOTES:

- Additional terms and conditions may apply. See GF's General Terms and Conditions on its website: www.grahamfield.com. Freight claims must be notated on the appropriate shipping documents and must be made with immediacy. International, federal and state regulations govern specific requirements for freight claims. Failure to abide by those regulations may result in a denial of the freight claim. GF will assist you in filing the freight claim. 2)
- Claims for any short shipment must be made within three (3) days of the invoice date.

