



Formulated with Chlorine to deodorize blood and body fluid spills.

COMBAT SPILLS

Control and Solidify Blood and Body Fluid Spills with Ease.



Red Z[®] Spill Control Solidifier

First developed by Safetec in the late 1980s, Red Z[®] Solidifier was used to help prevent infection during the AIDS epidemic. Unexpected spills can happen anytime, anywhere – no matter what type of industry. With the use of Red Z[®], you can reduce the risks from improper spill cleanup by relying on the solidifying agent to do the hard work for you.



Fast-Acting



Deodorizes



No Splashing

- Fast-acting and easy-to-use
- Eliminates splashing by solidifying and turning spills into a semi-solid (gel) mass
- Deodorizes spills (contains chlorine)
- Cost-effective cleanup solution

1oz (approx. 30cc) : **1L** (1,000cc)*



*1oz solidifies 1L of water-based fluid

Perfect for Unexpected Spills in:

- Healthcare & Laboratories
- Occupational Health & Safety
- Correctional Facilities
- Industrial & Manufacturing
- Transportation
- Schools
- Restaurants



Safetec of America, Inc.
1-800-456-7077 | safetec.com

Directions for Use

1. Wear appropriate protective clothing.
2. Sprinkle Red Z[®] evenly over spill. Allow sufficient time for solidification to occur.
3. Remove solidified material with scoop and scraper (or similar tool) and place in an appropriate container*.
4. Clean and disinfect contaminated areas after use and wash thoroughly.



Wear Protective Clothing



Sprinkle Red Z[®] Over Spill



Remove Solidified Spill

Expiration Dating

Red Z[®] is shelf stable and does not require an expiration date. In order for the product to remain shelf stable, it must be kept in a cool, dry place and in its original packaging. If contents are hard or compact, simply shake the product to loosen, before use.

*Dispose of infectious and hazardous waste in accordance with all local, state and federal regulations.

Red Z[®] Order Information




Order No.	Description	Case Qty.
41101	Red Z [®] 5 oz. shaker top bottle	24/case
41099	Red Z [®] 7.75 oz. shaker top bottle	24/case
41102	Red Z [®] 11 oz. shaker top bottle	12/case
41103	Red Z [®] 15 oz. shaker top bottle	12/case
41109	Red Z [®] .75 oz. single-use pouch	288/case
41117	Red Z [®] 10g single-use pouch	200/case
41119	Red Z [®] 21g single-use pouch	100/case
41120	Red Z [®] 2 oz. single-use pouch	100/case
41115	Red Z [®] 3.5lb. bucket	1 each
41116	Red Z [®] 17.5lb. bucket	1 each
41220	Red Z [®] 50lb. bucket	1 each

Kits Containing Red Z[®]



Order No.	Description	Case Qty.
17120	EZ-Cleans Spill Kit (poly bag refill)	24 kits
17121	EZ-Cleans Plus [®] Kit (poly bag refill)	24 kits
15201	Multi-Purpose Spill Kit (plastic case)	1 kit
25000	National Standard EZ-Cleans Kit (plastic case)	24 kits
25005	National Standard EZ-Cleans Kit (plastic case)	12 kits
25006	National Standard EZ-Cleans Kit (metal case)	24 kits
17100	UP Compliance Kit (poly bag refill)	24 kits
17102	UP Compliance Kit (plastic case)	12 kits
15205	Spill Leader Kit (plastic case)	1 kit



OSHA Compliance 1910.1030 – Bloodborne Pathogens Standard

Universal Precautions is an approach to infection control. According to the concept of Universal Precautions, all human blood and certain human body fluids are treated as if known to be infectious.