
meeting all the requirements of U.S. OSHA 29 CFR 1910.134 such as medical evaluation, training and fit testing. Select and use respirator in accordance with all applicable regulations, standards and professional guidance. Fit testing must be performed while the test subject is wearing any applicable safety equipment that may be worn during actual respirator use which could interfere with respirator fit, such as hair bonnets and eyewear. In Canada, CSA standard Z94.4 requirements must be met and/or requirement of the applicable jurisdiction, as appropriate. Follow all applicable local regulations. This respirator is designed for occupational/professional use by adults who are properly trained in its use and limitations. The 3M FT-10 (sweet solution) or FT-30 (bitter solution) Qualitative Fit Test Apparatus or other OSHA accepted fit test protocols are recommended for fit testing this respirator.

- Inspect respirator before each use to ensure that it is in good operating condition. Examine all the respirator parts for signs of damage including the two headbands, staples, nose foam, nose clip and facepiece material. The respirator should be disposed of immediately upon observation of damaged or missing parts.
- Leave the contaminated area immediately and contact your supervisor if dizziness, irritation or other distress occurs.
- Respirator may be used until damaged, breathing becomes difficult or contaminated with blood or body fluids. Discard after every use when used for surgical procedures. Follow national, state, local, and facility infection control guidance and policies.

Fitting Instructions: Must be followed each time respirator is worn.

- Remove respirator from packaging and hold with straps facing upward (Fig. 1a). Place the bottom strap under the center flaps next to the "ATTENTION" statement (Fig. 1b).
- Fully open top and bottom panels, bending nosepiece around thumb at center of foam. Straps should separate when panels are opened. Make certain bottom panel is unfolded and completely opened (Fig. 2).
- Place respirator on your face so that the foam rests on your nose and the bottom panel is open under chin. Pull the top strap over your head and position it high on the back of the head (Fig. 3a). Then pull the bottom strap over your head and position it around the neck and below the ears (Fig. 3b). Adjust for a comfortable fit by pulling top panel toward the bridge of the nose and bottom panel under chin.
- Place your fingertips from both hands at the top of the metal nosepiece. Using two hands, mold the nose area to the shape of your nose by pushing inward while moving your fingertips down both sides of the nosepiece (Fig. 4).
- ⚠ Pinching the nosepiece using one hand may result in improper fit and less effective respirator performance. (Use two hands.)**
- Perform a User Seal Check (Fig. 5). To check the respirator-to-face seal, cover the middle panel with one or both hands. Inhale and exhale sharply. Be careful not to disturb the position of the respirator. If air leaks around the nose, re-adjust the nosepiece as described in Step 4. If air leaks around respirator edges, adjust panels and position of straps and make certain respirator edges fit snugly against the face. **If you CANNOT achieve a proper seal, DO NOT enter the contaminated area. See your supervisor.**

Removal Instructions:

- Without touching the respirator facepiece, slowly lift the bottom strap from around your neck up and over your head (Fig. 1).
- Then lift off the top strap (Fig. 2).
- Store or discard according to your facility's infection control policy (Fig. 3). Dispose of used product in accordance with applicable regulations.

Storage Conditions and Shelf Life:

Before use, store respirators in the original packaging, away from contaminated areas, dust, sunlight, extreme temperatures, excessive moisture and damaging chemicals. When stored in original packaging between temperatures from -4°F (-20°C) to +86°F (+30°C) and not exceeding 80% RH, the product may be used until the date specified on packaging located next to the "Use by Date" symbol.

Explanation of Symbols

	= Use by Date
	= Manufacturer
	= Date of Manufacture
	= Batch code
	= Consult Instructions for use

NIOSH Approved N95 Respirator:

At least 95% filtration efficiency against solid and liquid aerosols that do not contain oil.

FOR MORE INFORMATION

and assistance on 3M Health Care products, contact your local 3M representative or call 3M toll free in USA at 1-800-228-3957. In Canada, call 3M Helpline at 1-800-267-4414.

IMPORTANT:

Before use wearer must read and understand these User Instructions. Keep these User Instructions for reference.

⚠ WARNING

This respirator helps protect against certain particulate contaminants but does not eliminate exposure to or the risk of contracting any disease or infection. Misuse may result in sickness or death. For correct use, consult your supervisor and the User Instructions, or call 3M at 1-800-228-3957. In Canada, call 3M Helpline at 1-800-267-4414.

IMPORTANT: Before use wearer must read and understand these User Instructions. Keep these User Instructions for reference.

98-0060-0211-1_3
34-8727-3783-7
3100030292

Description:

The 3M™ Aura™ Health Care Particulate Respirator and Surgical Mask 1870+ is designed to help provide respiratory protection for the wearer. This respirator has a filter efficiency level of at least 95% against particulate aerosols free of oil[†]. It is fluid resistant[‡] and meets > 99% bacterial filtration efficiency (BFE)[§]. It is disposable and is cleared to be worn in surgery. It can fit a wide range of face sizes.

Not made with natural rubber latex.

Intended Use:

This product meets CDC guidelines for *M. tuberculosis* exposure control. As a respirator, it is intended to help reduce wearer exposure to certain airborne particles, including those generated by electrocautery, laser surgery, and other powered medical instruments. As a surgical mask, it is designed to be fluid resistant to splash and spatter of blood and other infectious materials[¶].

Contraindications:

3M recommends that this respirator is not for industrial use. Not for use with beards or other facial hair or conditions that prevent a good seal between the face and the sealing edge of the respirator. Does not protect against gases or vapors (i.e. anesthetic gases such as isoflurane or vapors from sterilants such as glutaraldehyde). **This respirator was not designed to be used by children.**

[†] Tested against 0.3 micron particles (mass median aerodynamic diameter) per U.S. 42 CFR 84.

[‡] Meets ASTM Fluid Resistant Challenge F1862 at 160mm Hg

[§] Meets ASTM Standard Test Method for Bacterial Filtration Efficiency F2101.

Use Instructions:
1. In the U.S., before occupational use of this respirator, a written respiratory protection program must be implemented

3M™ Health Care Particulate Respirator and Surgical Mask 1870+

User Instructions

IMPORTANT: Keep these User Instructions for reference

⚠ WARNING

This respirator helps protect against certain particulate contaminants but does not eliminate exposure to or the risk of contracting any disease or infection. Misuse may result in sickness or death. For correct use, consult your supervisor and the User Instructions, or call 3M at 1-800-228-3957. In Canada, call 3M Helpline at 1-800-267-4414.

IMPORTANT: Before use wearer must read and understand these User Instructions. Keep these User Instructions for reference.

98-0060-0211-1_3
34-8727-3783-7
3100030292

Description:

The 3M™ Aura™ Health Care Particulate Respirator and Surgical Mask 1870+ is designed to help provide respiratory protection for the wearer. This respirator has a filter efficiency level of at least 95% against particulate aerosols free of oil[†]. It is fluid resistant[‡] and meets > 99% bacterial filtration efficiency (BFE)[§]. It is disposable and is cleared to be worn in surgery. It can fit a wide range of face sizes.

Not made with natural rubber latex.

Intended Use:

This product meets CDC guidelines for *M. tuberculosis* exposure control. As a respirator, it is intended to help reduce wearer exposure to certain airborne particles, including those generated by electrocautery, laser surgery, and other powered medical instruments. As a surgical mask, it is designed to be fluid resistant to splash and spatter of blood and other infectious materials[¶].

Contraindications:

3M recommends that this respirator is not for industrial use. Not for use with beards or other facial hair or conditions that prevent a good seal between the face and the sealing edge of the respirator. Does not protect against gases or vapors (i.e. anesthetic gases such as isoflurane or vapors from sterilants such as glutaraldehyde). **This respirator was not designed to be used by children.**

[†] Tested against 0.3 micron particles (mass median aerodynamic diameter) per U.S. 42 CFR 84.

[‡] Meets ASTM Fluid Resistant Challenge F1862 at 160mm Hg

[§] Meets ASTM Standard Test Method for Bacterial Filtration Efficiency F2101.

Use Instructions:
1. In the U.S., before occupational use of this respirator, a written respiratory protection program must be implemented

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Exempt of latex of caoutchouc naturel.

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© 3M 2021