

Grasping Cuff

#9464

Figure A

Position cuff on back of hand, centering proximal interphalangeal joints (PIP's) under cup, hook and loop pads at back of hand, metal loop on palmar side.

Wrap around wrist, secure with hook and loop.




Figure B

To apply palmar piece, open hook and loop Dycem wraps, position on hand, palm up. See illustration. Note that metal loop is positioned by thumb when on left hand, and be the little finger when on the right.

Figure C and D

Place handle to be held at preferred grab point. Noting grab point, remove palmar piece from hand and wrap the hook and loop Dycem wraps tightly around the handle (See Figure D).

Figure E

Return handle with palmar piece attached to palm of hand. Attach hook and loop strap from palmar piece to hook and loop pad at back of cuff. Pass hook and loop strap through metal loop on palmar piece, double back and also fasten to back of cuff. Insert cinching strap through large metal loop on wrist. Secure at desired tightness.

Figure F

For occasions when full finger wrap is not necessary, insert alternate strap through rectangular metal loop on palmar piece. Attach hook on palmar piece. Attach hook on strap to loop on palmar piece.

Fold strap back, adjust for proper snugness and secure.

