Introcan Safety® 3 Closed IV Catheter


en IV catheter designed to minimize inadvertent needlesticks, made of PUR, radiopaque, with fixation wings.

 For single use only	 Do not resterilize	 Sterilized using Ethylene Oxide.	 Non-pyrogenic	 Consult Electronic Instructions for Use. www.bbraunusa.com
 Use-by Date	 Reorder number	 Lot number	 Single Sterile Barrier System	 Single Sterile Barrier System with Protective Packaging Outside


Contents sterile in unopened, undamaged package. Not made with natural rubber latex, DEHP or PVC. Pictorials are for reference only. Pat.: www.bbraunusa.com/patents

RADIOPAQUE

Rx only P-8709-3 REV. 2/22

15328249

en

Fig. 1: Before needle removal



Fig. 2: After needle removal

en Description

The Introcan Safety® 3 Closed IV Catheter is designed to:

- Minimize inadvertent needlestick injuries with a passive safety needle-shielding feature.
- Improve catheter stability with an integrated stabilization platform designed to minimize movement within the vessel.
- Control the flow of blood which aids in the prevention of blood exposure.

The catheter may be used for any patient population with consideration given to adequacy of vascular anatomy and appropriateness for the solution being infused and duration of therapy.

Materials used

ABS, Chrome-Nickel Steel, POM, PP, PUR, Elastomer

Indications

- The Introcan Safety® 3 Closed Intravascular Catheter is inserted into a patient's vascular system for short term use to sample blood, monitor blood pressure, or administer fluids and blood intravascularly.
- The catheters may be used with power injectors at a maximum pressure of 325 psi with a luer lock connection only.

Contraindications

Introcan Safety® 3 Closed IV Catheter should not be used in patients with known hypersensitivity to any of the materials employed.

Guidelines

- This device features a multiple blood control septum. After needle withdrawal, blood from catheter hub is restricted until connection of a luer access device which opens the fluid path again.
- Observe Standard Precautions on ALL patients. Aseptic technique, proper skin preparation, and continued protection of the insertion site are essential.

- Use of only ISO Luer-slip and Luer-lock connectors ensures proper functioning.
- This device is intended for single use only and is provided sterile and non-pyrogenic.
- This device is sterile unless package has been opened or damaged.

Risks & Warnings

- Re-use of single-use devices poses a potential risk to patient and user. It may lead to loss of function of the device, or contamination which could cause severe injury or death.
- Always prime the catheter prior to arterial puncture.
- Verify adequate collateral circulation prior to arterial puncture.
- Non-ischemic complications of arterial catheterization, or even repeated attempts at catheterization, including bleeding, pseudoaneurysms, arteriovenous fistulae, nerve palsy, infection and injury to the tendon sheaths and adjacent structures during insertion.
- Arterial catheterization risks include injury to the adjacent structures during insertion, such as nerve and tendon injury, infection, ischemic injury such as vascular spasm and thromboembolism.
- Always clearly mark arterial lines to avoid accidental intra-arterial injection.
- Care must be taken to avoid needlesticks.
- Do not attempt to override or defeat the safety shield.
- In the unlikely event that the safety mechanism does not engage, keep the needle point away from the body and fingers at all times and immediately dispose into an approved sharps container.
- In the case of an unsuccessful IV start, remove the needle first to activate safety mechanism, then remove catheter from patient and discard both.
- Report needlesticks immediately and follow established institutional protocols.
- Do not bend the catheter/needle during insertion, advancement, or removal of the needle.
- Extreme care should be taken not to cut the catheter and possibly cause an embolism.
- Do not use scissors or sharp instruments at or near the insertion site.
- Never reinsert the needle inside the catheter once the needle has been partially or completely withdrawn as it may pierce and/or sever the catheter.
- Luer-slip connections should not be left unattended due to potential of disconnection.
- If blood is present inside the catheter, flush device according to institutional protocol.
- Not intended to facilitate or aid in the placement of vascular access devices such as guidewires, indwelling central venous catheters (CVC), peripherally inserted central catheters (PICC) and midline catheters into the vascular system.

Pressure Injection Information/Warnings and Cautions

- Prior to use with power injectors, ensure that a secure connection exists between the catheter and power injector.
- Patency of the catheter must be assured immediately before power injecting.
- Measures should be taken to avoid kinking or obstructing the catheter system during power injection to avoid product failure.
- Pressure limit settings on high pressure injector equipment may not prevent over pressurizing an occlude or partially occluded catheter, which may cause catheter failure.
- Follow the contrast media manufactures instructions for use, contraindications, warnings and precaution.
- Due to variations in add-on devices, tubing, injectate, temperature and pressure limit setting the listed injection flow rates may not be achievable.

Power Injection Flow Rates

Gauge / Length	Flow Rate (mL/s) Contrast Media Viscosity at 20°C 2.3 mPa*s	Flow Rate (mL/s) Contrast Media Viscosity at 20°C 27.5 mPa*s
24 Ga. x 3/4 in.	5.0	2.5
22 Ga. x 1 in.	8.0	3.5
20 Ga. x 2 in.	10.5	4.0
20 Ga. x 1 in.	10.5	4.0
20 Ga. x 1 - 1/4 in.	10.5	4.0
18 Ga. x 1 - 1/4 in.	14.0	5.0
18 Ga. x 1 - 3/4 in.	14.0	5.0
16 Ga. x 1 1/4 in.	15.5	5.5
16 Ga. x 2 in.	15.5	5.5
14 Ga. x 1 1/4 in.	16.5	5.5
14 Ga. x 2 in.	16.5	5.5

Duration of Use

- Change device according to CDC Guidelines and/or institutional protocol.
- The insertion site should be checked at regular intervals.
- Introcan Safety® 3 Closed IV Catheter should be removed in the event of local or systemic signs of infection.

Instructions for Use: Use Aseptic Technique. Do not rotate catheter hub prior to insertion.

1. Gather the necessary supplies.
2. Select and prepare site per institutional protocol.
3. Apply tourniquet.
4. Holding flashback chamber, remove protective guard in a straight outward motion. Inspect device.
5. Hold the device by the flashback chamber with thumb and fingers on opposite sides. Verify push-off plate and needle bevel are in the "up" position.
6. Confirm catheter hub is properly seated on flashback chamber.
7. Anchor the vessel with gentle skin traction. Adjust angle of insertion and puncture a suitable vessel. If vessel puncture was successful, blood will immediately be visible inside the transparent flashback chamber (see Fig. A).
8. Lower and advance the entire catheter and needle unit slightly (approximately 1/8 in. or 3 mm) to ensure catheter tip is in the vessel (see Fig. B).
9. Using the push-off plate, advance the catheter off the needle and observe for blood return between the needle and catheter (see figure C). After confirmation, continue advancing catheter off the needle and into the vessel.
10. Remove tourniquet.
11. Gently press the stabilization platform to the skin to stabilize the catheter (see Fig. D).
12. Stabilize the catheter hub to prevent catheter dislodgement during needle removal. Withdraw the needle straight out with a controlled and continuous motion (minimize rotation or bending of the needle). The metal safety shield will automatically attach to needle tip as needle tip exits catheter hub. Blood flow from catheter hub is restricted after needle is removed and until connection of a Luer access device is made (see Fig. E).
13. Dispose of the shielded needle immediately into an approved sharps container.
14. Immediately connect the infusion line or accessory device and cover puncture site with a sterile and transparent dressing per institutional protocol (see Fig. F).
15. After removal, dispose of the catheter per institutional protocol.

Distributed by:
B. Braun Medical Inc.
Bethlehem, PA 18018-3524 USA
1-800-227-2862

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B. Braun Melsungen AG