Rolyan® Progressive® Elbow/Knee Splint:

A812-410 through A812-414

Figures referred to in the instructions are on page 6.

INDICATIONS

- For mild to moderate spasticity or contractures
- Can be used in the presence of arthritis, dementia, Parkinson's disease, and osteoporosis

CONTRAINDICATIONS

- Not for use with severe fixed deformities
- Not for use with severe edema

PRECAUTIONS

- **This splint is to be fitted initially by a healthcare professional who is familiar with the purpose for which it is prescribed. The healthcare professional is responsible for providing wearing instructions and precautions to other healthcare practitioners, care providers involved in the patient's care, and the patient.**
- **If unusual swelling, skin discoloration or discomfort occurs, use should be discontinued and a healthcare professional consulted.**
- **Be careful that the splint is not so tight that circulation is impaired.**
- **After heating splinting material, always check its temperature before applying it to the patient.**

CARE AND CLEANING

- Splint should be kept away from open flames.
- Stay should be removed and all hook-and-loop closures fastened before hand washing and air drying the soft splint.
- Be sure the splint pieces are completely dry before reapplying them.
- The stay will lose its shape in temperatures over 135°F (57°C).

INSTRUCTIONS FOR USE

1. Place splint onto extremity to check for fit. Splint should allow unobstructed motion of shoulder or hip and clear bony protrusions. When wrapped around extremity, edges of sections should overlap approximately 2" (5 cm). Excess can be trimmed to reduce bulk.
2. Without allowing them to touch each other, heat the three parts of the stay in water heated to 160° to 170°F (70–75°C) until soft enough to mold. Remove from heat pan and place on a towel on a flat surface.

Forming the stay

Note that the **Ezeform®** stay pieces will bond quickly, so place them carefully.

3. Place largest piece on flat surface.
4. Line up top edge of cross bar of T-shaped piece with top edge of cross bar of largest piece (**A**). Cutout area of T-shaped piece should line up with extended top of largest piece. Press pieces together gently.

5. Place straight piece so top is even with cutout edge of T-shaped piece and centers line up (**B**).
 6. Bend extended top of largest piece over other two pieces (**C**) to fit into cutout area of T-shaped piece.
 7. Fold two slit sides of largest piece onto middle of stay, forming a "T" (**D**).
 8. If necessary, trim edges of stay so they line up evenly.
 9. If necessary, reheat stay. Place stay on extremity with center seam away from skin and cross bar on upper arm or thigh. Mold to appropriate degree of elbow or knee extension (**E**).
 10. For greater comfort, flare ends of stay slightly.
- ### Applying the splint
11. Open fastener on stay pocket and put stay into pocket (**F**). Close flap over cross bar of stay. Top D-ring strap also has hook-and-loop closure to keep it in place on pocket flap.
 12. Place splint on extremity and secure D-ring straps. A properly applied splint for the upper extremity is shown in **G** and **H**.
 13. As condition changes, splint may be progressed by heating stay and remolding it.

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Patterson Medical Ltd.