

# Rolyan® Incremental Wrist Hinge:


A352-100 and A352-101

**A**

**B**

**C**

**D**

A pattern for forearm and hand splint bases that can be used with the Rolyan Incremental Wrist Hinge is included.

## INDICATIONS

- To statically position wrist or block out undesired range of motion

## CONTRAINDICATIONS

- Not in the presence of unstable wrist fractures
- Not in the presence of fixed bony contractures
- Not for moderate to severe flexor tone

## PRECAUTIONS


- This product is to be fitted initially by a healthcare professional who is familiar with the purpose for which it is prescribed. The healthcare professional is responsible for providing wearing instructions and precautions to other healthcare practitioners, care providers involved in the patient's care, and the patient.**
- If unusual swelling, skin discoloration or discomfort occurs, use should be discontinued and a healthcare professional consulted.**
- Be careful that the straps are not so tight that they interfere with circulation.**
- After heating splinting material, always check its temperature before applying it to the patient.**

## CARE AND CLEANING

Hinge can be cleaned with a damp cloth.

## INSTRUCTIONS FOR USE

Splint components are referred to by the following names:

## Forming splint bases

- Fabricate two-part splint, which must consist of the following:
  - forearm base that extends proximal to ulnar styloid on dorsal surface of forearm

- hand base (circumferential for most applications)
- Axis of hinge must be aligned just distal to peak of radial styloid, and both splint portions must be large enough and extend far enough radially to accommodate hinge (A).  
The Rolyan Tenodesis Flexor Tendon Repair Pre-Cut Splint Blank is recommended for use following repair of a flexor tendon injury.

## Attaching wrist hinge to splint bases

- Apply splint bases.
- With contoured leg of hinge distal, position hinge along radial border in alignment with anatomical axis of wrist (just distal to peak of radial styloid).
- Trace slots in legs of hinge onto both splint bases.
- Remove from patient and use a  $\frac{3}{16}$ " drill bit or punch to make four holes within marked areas, two in hand base and two in forearm base (B). The holes must be at least  $\frac{1}{2}$ " (1.3 cm) apart to accommodate the 2-part screws without overlapping.
- Insert post sections of 2-part screws, with flat surface inside splint, into holes. Press until heads are flush with splint.
- Place wrist hinge onto splint over post sections. From outside of hinge, screw slotted-head sections through slots in legs of hinge and into post sections (C). Tighten with screwdriver. Slip-resistant material can be used to help hold post sections while tightening.

## Joining ulnar connecting pieces

- Mark position for ulnar hole in hand base. (See splint pattern on page 2 for proper hole position.) Using a  $\frac{3}{16}$ " drill or punch, make hole.
- If base is padded, cut away padding around hole.
- With splint on patient, position ulnar connecting piece of hand base over ulnar connecting piece of forearm base. Be sure that connecting pieces line up with axis of wrist.
- Mark connecting piece of forearm base through hole in hand base. Make sure that patient will be able to flex and extend wrist when connecting pieces are attached at this point.
- Remove splint and use drill or punch to make two holes that overlap slightly, forming a figure 8.
- Put post section of 2-part screw through forearm base, with flat head inside splint.


- Place hole in hand base over post section in connecting piece of forearm base. Put slotted-head screw into post section and tighten with screwdriver.

## Adjusting range of motion

Note that degrees marked on hinge represent degree at which hinge is set, not necessarily degree at which wrist is positioned.

### 15. For a static hinge

- Line up arrow on hinge with desired amount of flexion or extension. Place adjustment screw through that hole, making sure that it goes through window in back of hinge.
- As the patient gains flexion or extension, simply loosen adjustment screw and

progress to next degree level.

### 16. For free-moving hinge with end-range blocking

- Line up arrow on hinge with greatest amount of extension desired. Place adjustment screw in hole that is marked 20° greater than hole to which arrow points (D). For example, if hinge must stop at 30° of extension, point arrow at 30°, and place adjustment screw at 50°.

**Optional splint pattern**  
**Optioneel patroon voor spalk**  
**Modèle d'attelle en option**  
**Optionales Schienenmuster**  
**Modello per splint opzionale**  
**Patrón opcional de la férula**


ulnar connecting pieces  
 ulnaire verbindingsstukken  
 pièces de connexion cubitales  
 ulnare Verbindungsstücke  
 porzioni di collegamento ulnare  
 piezas de conexión ulnares

mark circle onto  
 splint blank (a hole is  
 punched here)  
 cirkel aangeven op de  
 spalkdeken (hier wordt  
 het gat gemaakt)  
 tracer un cercle sur  
 l'ébauche d'attelle (un  
 trou y sera percé)  
 Kreis auf  
 Schienenplatte  
 aufzeichnen (dort wird  
 ein Loch gestanzt)  
 contrassegnare il  
 cerchietto sullo splint  
 da modellare (qui  
 viene praticato un foro)  
 marque un círculo en  
 el patrón para férula  
 (aquí se perfora un  
 orificio)


**Hand Base**  
**Spalkdeken**  
**Base de main**  
**Hand-Grundkörper**  
**Base per la mano**  
**Base de la mano**

palmar section of splint  
 palmair deel van de spalk  
 section palmaire de  
 l'attelle  
 palmarer Teil der Schiene  
 sezione palmare dello  
 splint  
 sección palmar de la  
 férula

**Forearm Base**  
**Onderlegger onderarm**  
**Base d'avant-bras**  
**Unterarm-Grundkörper**  
**Base per l'avambraccio**  
**Base del antebrazo**

**For left hand**  
**Voor de linkerhand**  
**Pour la main gauche**  
**Für die linke Hand**  
**Per mano sinistra**  
**Para mano izquierda**

**For right hand**  
**Voor de rechterhand**  
**Pour la main droite**  
**Für die rechte Hand**  
**Per mano destra**  
**Para mano derecha**

## Optional splint pattern for use with Rolyan Incremental Wrist Hinge

Pattern is on page 2.

### INSTRUCTIONS FOR USE

The pattern fits small- to medium-sized people. It can be reduced or enlarged on a photocopier to accommodate a variety of hand sizes.

1. Photocopy pattern and cut it out. Place both pieces onto patient to check for proper fit.
2. Trace pattern onto splinting material. Heat splint blanks until soft and cut them out.
3. If desired, pre-pad the splint bases with 1/8" (3.2mm) padding material.

### Molding forearm splint base


4. Heat forearm base until moldable.
5. Mold on dorsal surface of forearm so ulnar connecting piece crosses axis of wrist and ends at distal wrist crease.
6. Check for fit and trim as needed.
  - Lengthwise, splint should extend approximately two-thirds of length of forearm.
  - Widthwise, splint should extend only to radial and ulnar mid-lines of forearm—except ulnar connecting piece, which may extend past ulnar mid-line.

- Distal edge of ulnar connecting piece should end at distal wrist crease.
- Ulnar connecting piece and splint body should clear ulnar styloid (note that connecting piece must be at least 1" [2.5 cm] wide to accommodate connecting screw).

### Molding hand splint base

7. Apply desired straps.
8. Heat hand base until moldable.
9. With forearm splint base on patient and wrist in as much flexion as possible, position hand base with its ulnar connecting piece overlapping connecting piece on forearm base by at least 1" (2.5 cm) (refer to diagrams on page 2 for proper orientation of splint bases).
10. Form hand base through thumb web space.
11. Check for fit and trim as needed.
  - Palmar section should end on volar surface mid-way across palm.
  - Palmar section should clear distal palmar crease for unimpeded MCP flexion.

- Splint should allow desired thumb range of motion.
  - Splint must be long enough to accommodate hinge.
12. Attach desired straps and wrist hinge.
  13. Attach 2-part connecting screw to ulnar side of splint bases.

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Patterson Medical, Ltd.**