**CRAMER ATHLETIC SPORTS TAPE**

**Cramer 750, Cramer 950, Cramer 1050i, Cramer Non-Tear Stretch Tape, Cramer AT Pro**

**Instructions For Use**

- Wash the area to be taped with soap and water and towel dry.
- (Optional) Apply Cramer Tape Underwrap over the skin, being careful to keep underwrap smooth. Use Cramer Heel & Lace pads over blister-prone areas if needed.
- Apply Cramer Athletic Tape to the area. Wrap it securely, keeping the tape smooth, avoiding wrinkles.
- Remove tape immediately after use to avoid skin irritation. Use bandage scissors or tape cutters to cut tape.
- Wash area with soap and water. Caution: Skin irritation or blisters may occur on area of application due to friction while using and removing tape.

**Cramer Pro-Lastic Tear Stretch Tape**

**Instructions For Use**

- If used directly on the skin, wash the area to be taped with soap and water and towel dry.
- Apply Cramer Pro-Lastic Tape to the area. Wrap it securely, keeping the tape smooth, avoiding wrinkles.
- Remove tape immediately after use to avoid skin irritation. Use bandage scissors or tape cutters to cut tape.
- Wash area with soap and water if used directly on skin.
- Caution: Skin irritation or blisters may occur on area of application due to friction while using and removing tape.

**Cramer Eco-Flex Tape**

**Instructions For Use**

- If used directly on the skin, wash the area to be taped with soap and water and towel dry.
- Apply Cramer Eco-Flex Tape to the area. Wrap it securely, keeping the tape smooth, avoiding wrinkles. Cut tape to desired length.
- Remove tape immediately after use to avoid skin irritation. Use bandage scissors or tape cutters to cut tape.
- Caution: Skin irritation or blisters may occur on area of application due to friction while using and removing tape.

**Cramer Tape Underwrap**

**Instructions For Use**

- Use Cramer Tape Underwrap in conjunction with Cramer Athletic Tape.
- Wash the area to be taped with soap and water and towel dry.
- Apply Cramer Tape Underwrap over the skin, being careful to keep underwrap smooth.

**Intended Use: Athletic tape assists in the prevention and healing of sprains, strains and contusions by providing support and limiting movement of muscles and joints.**

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