
HemoCue Hb 201+
Operating Manual
Manuel d'utilisation

HemoCue Hb 201+ *Operating Manual*

HemoCue Hb 201+ *Manuel d'utilisation*

US

The HemoCue Hb 201+ analyzer provides determination of hemoglobin quickly, easily and with lab quality results. Capillary, venous or arterial whole blood may be used. This manual provides the basic instructions for use as well as technical specifications. Additional information may be obtained from HemoCue, Inc.

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Components Composants

US

1. The HemoCue Hb 201+ analyzer
2. AC adapter
3. 4 type AA batteries*
4. A vial of HemoCue Hb 201 Microcuvettes*
5. Individually packaged HemoCue Hb 201 Microcuvettes*

The HemoCue Hb 201+ and its adapter are delivered in a carton. Open the carton on a stable surface and lift out the instrument and accessories. On the bottom of the analyzer there is a lid covering the battery compartment. Press the flap to remove the lid. Place 4 type AA batteries in the battery compartment, observing the indication of polarity in the battery holders. Replace the lid. Consult local environmental authorities for proper disposal of batteries.

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
Start-up Démarrage

US

- 1a. If AC power is available, plug the supplied power adapter into the power inlet at the back of the analyzer.
- 1b. If no power is available, insert the 4 type AA batteries, 1.5 V, into the battery compartment.
- 1c. If a battery symbol appears on the display, the batteries are running low on power. The analyzer will continue to give accurate results, but the batteries should be replaced as soon as possible.
2. Pull the cuvette holder out to the loading position.
3. Press and hold the left button until the display is activated (all symbols appear on the display).
4. The display shows the version number of the program, after which it will show "☼" and "Hb". During this time the analyzer will automatically verify the performance of the optronic unit by performing an automatic SELFTEST..
5. After 10 seconds, the display will show three flashing dashes and the HemoCue symbol. This indicates that the HemoCue Hb 201+ has passed the SELFTEST and is ready for use. If the SELFTEST fails, an error code will be displayed.

To turn the analyzer off, press and hold the left button until the display reads OFF and then goes blank. Disconnect the equipment from the power source by removing the AC adapter from the outlet.

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
Mesure Sang capillaire



US

1. To perform a test using capillary blood, the cuvette holder should be in its loading position. The display will show three flashing dashes and the HemoCue symbol.
2. Make sure the patient's hand is warm and relaxed. Use only the middle or ring finger for sampling. Avoid fingers with rings on.
3. Clean the finger with alcohol or a suitable disinfectant and allow to dry or wipe off with a dry, lint-free wipe.
4. Using your thumb, lightly press the finger from the top of the knuckle towards the tip. This stimulates the blood flow towards the sampling point.

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
US

5. For best blood flow and least pain, sample at the side of the fingertip, not in the center.
6. While applying light pressure towards the fingertip, puncture the finger using a lancet.
7. Wipe away the first 2 or 3 drops of blood.
8. Re-apply light pressure towards the fingertip until another drop of blood appears.

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

11


US

9. When the blood drop is large enough, fill the microcuvette in one continuous process. Do NOT refill!
10. Wipe off excess blood from the outside of the microcuvette with a clean, lint-free wipe, being careful not to touch the open end of the microcuvette, which could result in blood being drawn out of the microcuvette.
11. Look for air bubbles in the filled microcuvette. If present, discard the microcuvette and fill a new microcuvette from a new drop of blood. Small bubbles around the edge can be ignored.

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US

12. Place the filled microcuvette in the cuvette holder. This must be performed within ten minutes after filling the microcuvette!
13. Gently slide the cuvette holder to the measuring position.
14. During the measurement “” and three fixed dashes will be shown on the display.
15. After 15–60 seconds, the hemoglobin value of the sample is displayed. The result will remain on the display as long as the cuvette holder is in the measuring position. When operating on battery power, the analyzer will automatically turn off after approximately 5 minutes.
16. Although the reagents are present in the microcuvette in extremely low quantities, consult local environmental authorities for proper disposal. Always handle blood specimens with care, as they might be infectious.

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Measuring Control material, venous or arterial blood

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


④

US

The HemoCue Hb 201 system can be verified on the days of testing through the use of commercially available hematology controls. Contact HemoCue, Inc. for control information.

1. The blood or control material should be well mixed prior to performing a measurement.
- 2ab. Place a drop of blood or control material onto a hydrophobic surface, e.g., a plastic film, using a pipette or a suitable blood transfer device.
3. Fill the microcuvette in one continuous process. Do NOT refill! Wipe off excess blood from the outside of the microcuvette with a clean, lint-free wipe, being careful not to touch the open end of the microcuvette, which could result in blood being drawn out of the microcuvette.
4. Look for air bubbles in the filled microcuvette. If present, discard the microcuvette and fill a new microcuvette from a second drop of sample. Small bubbles around the edge can be ignored.

Perform the analysis as per steps 12–16 on page 15.

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Maintenance Entretien

US

The cuvette holder should be cleaned after each day of use.

1. Check that the analyzer is turned off. The display should be blank.
2. Pull the cuvette holder out to its loading position. Carefully press the small catch positioned in the upper right corner of the cuvette holder.
3. While pressing the catch, carefully rotate the cuvette holder towards the left as far as possible. Carefully pull the cuvette holder away from the analyzer.
4. Clean the cuvette holder with alcohol or mild detergent. It is important that the cuvette holder is completely dry before being replaced.

The optronic unit should be cleaned when directed to do so in the Troubleshooting Guide or as desired.

A dirty optronic unit may cause the analyzer to display an error code. To clean the optronic unit, proceed as follows.

5. Push the HemoCue Cleaner* swab into the opening of the cuvette holder. Move from side to side 5–10 times. If the swab is stained, repeat with a new swab. No further cleaning is required if the swab remains clean.
6. Wait 15 minutes before replacing the cuvette holder and using the analyzer.

The cover may be cleaned with alcohol or a mild soap solution.

*Not included. As an alternative to the HemoCue Cleaner, a cotton tip swab moistened with a alcohol (without additive) or water may also be used for cleaning.

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US Troubleshooting Guide

If you are unable to resolve the problem by following this Troubleshooting Guide, please contact HemoCue, Inc. The analyzer has no serviceable parts.

| Symptom | Explanation | Action |
|-----------------------------------|--|---|
| The analyzer shows an error code. | May be a temporary fault | Turn off the analyzer and turn it on again after 30 seconds. Take a new microcuvette and repeat the measurement. If the problem continues, see specific error code below. |
| E00 | No stable endpoint is found within the time range. 1. The cuvette is faulty. 2. The circuit board is out of order. | 1a. Check the expiration date for the microcuvettes. 1b. Take a new microcuvette and repeat the measurement. 2. The analyzer needs service. Contact HemoCue, Inc. |
| E01–E05 | 1. Dirty optronic unit or faulty electronic or optronic unit | 1a. Turn off the analyzer and clean the optronic unit as described in the maintenance section. 1b. The analyzer needs service. Contact HemoCue, Inc. |
| E06 | 1. Unstable blank value The analyzer might be cold. | 1. Turn off the analyzer and allow it to reach room temperature. If the problem continues, the analyzer needs service. Contact HemoCue, Inc. |
| E07 | 1. The battery power is too low. | 1a. The batteries need to be replaced. Turn off the analyzer and replace the batteries, 4 type AA. 1b. Use the power adapter. |
| E08 | The absorbance is too high. 1. An item is blocking the light in the cuvette holder. | 1a. Check that the analyzer and microcuvettes are being used according to the HemoCue Hb 201+ operating manual and instructions for use. 1b. The analyzer needs service. Contact HemoCue, Inc. |
| E09–E30 | 1. Dirty optronic unit or faulty electronic or optronic unit | 1a. Turn off the analyzer and clean the optronic unit as described in the maintenance section. 1b. The analyzer needs service. Contact HemoCue, Inc. |

| Symptom | Explanation | Action |
|---|--|--|
| HHH | 1. Measured value exceeds 25.6 g/dL (256 g/L, 15.9 mmol/L) | |
| No characters on the display | <ol style="list-style-type: none"> 1. The analyzer is not receiving power. 2. If on battery power, the batteries need to be replaced. 3. The display is out of order. | <ol style="list-style-type: none"> 1a. Check that the power adapter is connected to the AC power supply. 1b. Check that the power adapter is securely connected to the analyzer. 1c. Check that the cable is not damaged. 2. Turn off the analyzer and replace the batteries, 4 type AA. 3. The analyzer needs service. Contact HemoCue, Inc. |
| The display gives erroneous characters. | <ol style="list-style-type: none"> 1. The display is out of order. 2. The microprocessor is out of order. | <ol style="list-style-type: none"> 1. The analyzer needs service. Contact HemoCue, Inc. 2. The analyzer needs service. Contact HemoCue, Inc. |
| The display shows "▢▢". | <ol style="list-style-type: none"> 1. The batteries need to be replaced. 2. If on AC power, the power adapter or the circuit board is out of order. | <ol style="list-style-type: none"> 1. Turn off the analyzer and replace the batteries, 4 type AA. 2a. Check that the power adapter is properly connected and working. 2b. The analyzer needs service. Contact HemoCue, Inc. |
| The display does not switch from "⊠" and "Hb" to three flashing dashes and "🔍" (ready for measuring). | <ol style="list-style-type: none"> 1. The magnet in the cuvette holder may be missing. 2. The magnetic sensor is out of order. | <ol style="list-style-type: none"> 1. The analyzer needs service. Contact HemoCue, Inc. 2. The analyzer needs service. Contact HemoCue, Inc. |

| Symptom | Explanation | Action |
|---|--|---|
| <p>Measurements on control materials are out of range – either too HIGH or too LOW.</p> | <ol style="list-style-type: none"> 1. The microcuvettes are beyond their expiration date, damaged or have been improperly stored. 2. The optical eye of the microcuvette is contaminated. 3. The control has not been mixed properly and/or is not at room temperature. 4. Air bubbles in the microcuvette 5. The optronic unit is dirty. 6. The control is not suitable for use with the HemoCue Hb 201+ system. 7. The calibration of the analyzer has been changed. 8. The controls are beyond their expiration dates or have been improperly stored. | <ol style="list-style-type: none"> 1. Check the expiration date and the storage conditions of the microcuvettes. 2. Remeasure the sample with a new microcuvette. 3. Make sure that the control is mixed properly and at room temperature. 4. Check the microcuvette for air bubbles. Remeasure the sample with a new microcuvette. 5. Clean the optronic unit as described in the maintenance section. 6. Only use controls intended for the HemoCue Hb 201+ system. Contact HemoCue, Inc. for control information. 7. The analyzer needs service. Contact HemoCue, Inc. 8. Check the expiration date and the storage conditions of the control. Take a new microcuvette and repeat the measurement from a new vial/bottle of control. |
| <p>Measurements on patient samples are higher or lower than anticipated.</p> | <ol style="list-style-type: none"> 1. Improper sampling technique 2. The microcuvettes are beyond their expiration date, damaged or have been improperly stored. 3. The optical eye of the microcuvette is contaminated. 4. Air bubbles in the microcuvette 5. The optronic unit is dirty. 6. The calibration of the analyzer has changed. | <ol style="list-style-type: none"> 1. See pages 8-17 in this manual. 2. Check the expiration date and the storage conditions of the microcuvettes. Check the entire system with a commercial control. 3. Remeasure the sample with a new microcuvette. 4. Check the microcuvette for air bubbles. Remeasure the sample with a new microcuvette. 5. Clean the optronic unit as described in the maintenance section. 6. The analyzer needs service. Contact HemoCue, Inc. |

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Specifications

General

The HemoCue Hb 201+ is a system used for the determination of the total amount of hemoglobin in whole blood. The system consists of a specially designed analyzer with specially designed microcuvettes containing dried reagents. The microcuvette serves as pipette, reaction vessel and as a measuring microcuvette.

No dilution is required. The hemoglobin measurement takes place in the analyzer, which follows the progress of the reaction until the steady state has been reached. The system is factory calibrated against the hemoglobincyanide (HiCN) method, the international reference method for the determination of the hemoglobin concentration in blood.

Intended Use

Quantitative determination of hemoglobin in capillary, venous or arterial blood. The HemoCue Hb 201+ system is used for quantitative determination of hemoglobin in whole blood using a specially designed analyzer, the HemoCue Hb 201+ analyzer, and specially designed microcuvettes, the HemoCue Hb 201 Microcuvettes. HemoCue Hb 201 Microcuvettes are for In Vitro Diagnostic use only. The HemoCue Hb 201+ analyzer is only to be used with HemoCue Hb 201 Microcuvettes.

Theory

The reaction in the microcuvette is a modified azidemethemoglobin reaction. The erythrocyte membranes are disintegrated by sodium deoxycholate, releasing the hemoglobin. Sodium nitrite converts the hemoglobin iron from the ferrous to the ferric state to form methemoglobin, which then combines with azide to form azidemethemoglobin.

Reagents

40% w/w Sodium Deoxycholate, 18% w/w Sodium Azide, 20% w/w Sodium Nitrite, 22% w/w nonreactive ingredients.

Sample material

Capillary, venous or arterial blood may be used. Appropriate anticoagulants in solid form (eg EDTA, heparin or heparin/fluoride) may be used. Mix all samples thoroughly on a mechanical mixer for at least 2 minutes or invert the tube 8-10 times by hand. Alternatively, follow the local recommendations. Hemoglobin remains unchanged for days, provided that the blood does not become infected. If the specimen has been stored in a refrigerator, it will be viscous and the blood should be allowed to warm up to room temperature before mixing.

Storage and environmental requirements of the HemoCue Hb 201 Microcuvettes

Use the HemoCue Hb 201 Microcuvettes prior to their expiration date. The expiration date is printed on each package.

Storage of microcuvettes kept in a vial

The microcuvettes are to be stored at room temperature (59–86 °F, 15–30 °C). Do not refrigerate. The expiration date of the HemoCue Hb Microcuvettes in a sealed vial is printed on the vial. Once the seal is broken, the microcuvettes are stable for three months. Always keep the container properly closed.

Storage of the individually packaged cuvettes

The microcuvettes are to be stored at room temperature (59–86 °F, 15–30 °C). Do not refrigerate. The individually packaged microcuvettes are stable until the date printed on each package.

HemoCue Hb 201+ analyzer

The analyzer can be stored at 32–122 °F (0–50 °C). Operating temperature is 59–86 °F (15–30 °C). Allow the analyzer to reach ambient temperature before use. The analyzer should not be operated at high (i.e. > 90 % non-condensing) humidity.

Quality Control

The HemoCue Hb 201+ analyzer has an internal electronic "SELFTEST". Every time the analyzer is turned on, it will automatically verify the performance of the optronic unit of the analyzer. This test is performed every second hour if the analyzer remains switched on.

Upon passing the SELFTEST, the display will show the HemoCue symbol and three flashing dashes, indicating that the analyzer is ready to perform a measurement. An error code will be displayed if the SELFTEST fails. If use of liquid control material is required by local or other regulations, contact HemoCue, Inc. for control information.

Measuring range

0–25.6 g/dL (0–256 g/L, 0–15.9 mmol/L).

Results above 25.6 g/dL (256 g/L, 15.9 mmol/L) will be displayed as HHH. Values above 23.5 g/dL (235 g/L, 14.6 mmol/L) must be confirmed using a suitable laboratory method.

Limitations

HemoCue Hb 201 Microcuvettes are for in vitro diagnostic use only. The HemoCue Hb 201+ analyzer is only to be used together with HemoCue Hb 201 Microcuvettes. For further limitations of

the procedure, see the HemoCue Hb 201 Microcuvettes package insert.

Expected values

Adult Males 13.0–17.0 g/dL (130–170 g/L, 8.1–10.5 mmol/L)

Adult Females 12.0–15.0 g/dL (120–150 g/L, 7.4–9.3 mmol/L)

Infants, after neonatal period 11.0–14.0 g/dL (110–140 g/L, 6.8–8.7 mmol/L)

Children, two years to teenage: gradual increase to adult normals. Due to a wide range of conditions (dietary, geographical, etc) which affect normal values, it is recommended that each laboratory establish its own normal range.

Technical Specifications

Dimensions: 3.35 x 6.30 x 1.69 inches (85 × 160 × 43 mm)

Weight: 0.77 pounds (350 g) (with batteries installed)

4 batteries type AA

Power adapter: CE marked

Only use adapters recommended by HemoCue or equivalent, as listed on page 28.

Pollution degree: 2

Overvoltage category: II

Atmospheric pressure: 700 hPa to 1060 hPa.

Equipment not suitable for use in the presence of flammable mixtures.

The instrument is tested according to IEC 61010-1, First edition, 1990 with Amendment 1, 1992 and Amendment 2, 1995 (EN 61010-1: 1993 with A2: 1995), EN 60601-1-2 and complies with the IVD Medical Device Directive 98/79/EC.

Signal input and signal output

Accessory equipment connected to the analog and digital interfaces must be certified according to the respective UL standards (i.e. UL 60950-1 for data processing equipment and IEC 60601-1

for medical equipment). Furthermore all configurations shall comply with the system standard IEC 60601-1-1. Connections of additional equipment to the signal input or signal output connections is deemed to be "configuring a medical system", and therefore assumes responsibility for the systems compliance within the requirements of IEC 60601-1-1. If in doubt, contact HemoCue AB.

Warranty

The analyzer carries a 24-month warranty from the day of receipt. After the warranty period, service/repair is carried out at fixed prices.

Spare parts – Accessories

The following spare parts and accessories are available:

Power adaptor

Battery Lid

HemoCue Cleaners

Hard Carrying Case

Cuvette holder

HemoCue Lancets

HemoCue Hb 201+ Control Cuvette

Symbols

Attention, see instructions for use

CE mark

Refers to Type BF Applied Part, Per EN60601-1

Class II equipment

The Medical Equipment with respect to electrical shock, fire and mechanical hazards only in accordance with UL-2601-1 and related collateral IEC 60601-1 standards for which the products has been found to comply by UL.

Only valid within the European Community. Indicates separate collection for waste of electrical and electronic equipment.

References

1. Makarem, A. In clinical Chemistry: Principles and Technics, 2nd ed., Henry, R.J., Cannon, D.C., and Winkelman, J.W., Harper and Row, Hagerstown, M.D., 1974, pp. 1125–1147.
2. Dacie and Lewis, Practical Haematology Ninth edition.
3. Jacques Wallach, M.D., Interpretation of Diagnostic Tests.

Recommended Adapter or equivalent:

| Country | Type | Manufacturer | Input | Output |
|------------|----------------------|--------------|----------------|---------------|
| USA/Canada | Mod nr FE3515060D035 | FRIWO | 120 V AC/60 Hz | 6V DC, 350 mA |

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