# Rolyan® Preformed Dynamic Pronation/Supination Splint:

Elbow Component—A833-2 through A833-5

Wrist Component—A833-11, A833-12, A833-21, A833-22, A833-31 and A833-32

# Rolyan® Dynamic Pronation/Supination Kit:

A833-1

*Be sure to note which splint you are using and follow the instructions that apply to that splint. Unless otherwise stated, these instructions apply to both the preformed splint and the kit. Figures referred to in the instructions are on page 10.*

## INDICATIONS

- For increasing pronation or supination of the forearm
- May be used following wrist or elbow fracture, Colles' fracture, Smith's fracture, proximal radial head fracture, and radial-lunate and other partial fusions of the carpus

## CONTRAINDICATIONS

- Not before sufficient bone healing has occurred
- Not before pins have been removed
- Not with severe edema

## PRECAUTIONS

- **Use of this splint requires careful supervision by a healthcare professional. Adjustments in the tension should be made only by a healthcare professional.**
- **This splint is to be fitted initially by a healthcare professional who is familiar with the purpose for which it is prescribed. The healthcare professional is responsible for providing wearing instructions and precautions to other healthcare practitioners, care providers involved in the patient's care, and the patient.**
- **If unusual swelling, skin discoloration or discomfort occurs, use should be discontinued and a healthcare professional consulted.**
- **After heating splinting material, always check its temperature before applying it to the patient.**
- **Be careful that the straps are not so tight that they interfere with circulation.**
- **After making an adjustment in the amount of force provided, evaluate the splint to make sure the tension is appropriate.**

## CARE AND CLEANING

- The splint will lose its shape in temperatures over 135°F (57°C). It should be kept away from sources of heat such as ovens, hot water, open flames, and sunny car windows.
- The splint can be cleaned with soap and lukewarm water. The straps can be washed with soap and water, but the adhesive portion should not be moistened. The splint and straps should be allowed to dry thoroughly before reapplication of the splint.

## INSTRUCTIONS FOR USE

If using the kit, an elbow and a wrist splint must be fabricated (see instructions 1 through 16). If using the preformed elbow and wrist splint components, begin with instruction 17. The items included in the preformed splint and the kit are listed below.

### Rolyan Preformed Dynamic Pronation/Supination Splint—Elbow Component

The preformed elbow splint includes the following items:

- pre-padded splint of Polyform® splinting material

- two self-adhesive D-ring straps
- two strap pads
- one metal housing for the tubing
- two 2-part screws
- Allen key

### Rolyan Preformed Dynamic Pronation/Supination Splint—Wrist Component

The preformed wrist splint includes the following items:

- Polyform splint with attached metal housing for tubing
- foam splint liner
- one self-adhesive straight strap with hook piece
- two self-adhesive D-ring straps
- two strap pads

### Dynamic Pronation/Supination Kit

The kit includes the following items:

- tubing
  - two metal housings for the tubing
  - four 2-part screws
  - Allen key
- Additional items needed when using the kit:
- low-temperature splinting material for splints
  - closed-cell padding for splints or stockinette for wrapping the extremity
  - strapping materials
  - punch or drill with 1/4" bit
  - screwdriver
  - optional crescent wrenches or vise and dowel (see instruction 13)

## DYNAMIC PRONATION/SUPINATION KIT

### Fabricating the elbow splint

1. Measure and mark a splint that extends from two-thirds of the length of the humerus past the elbow to one-third the length of the forearm. It should be as wide as one-half of the circumference of the humerus and proximal forearm.
2. If desired, pad the splint blank. After heating it, form a posterior elbow splint with the elbow in 90° of flexion.

### Fabricating the wrist splint

3. Measure for the wrist splint as follows:
  - from distal palmar crease to one-half of the length of the forearm, allowing at least 1 1/2" (3.8 cm) between the elbow splint and the wrist splint
  - one-half of the largest circumference of the forearm
  - one-half of the circumference of the hand at the MCPs
4. If desired, cut padding to match but do not attach it to the splint blank.
5. On the splint blank and the padding, mark the thumb hole approximately 2" (5.1 cm) from the distal and radial edges.
6. Cut a hole in the padding large enough to accommodate the thumb and clear an area of the splint around the thumb that will be rolled back (see instructions 7 and 10).
7. Apply the padding to the splint blank. After heating it, cut the thumb hole in the splint blank and roll the edges of the opening to the outside surface of the splint.

8. If possible, place the patient's hand in supination to achieve a gravity-assisted position.
9. Stand in front of the patient and insert the thumb through the thumb opening. Extend the splint over the volar surface of the forearm. Mold it on the radial side of the thumb hole and continue to position the splint symmetrically on the forearm.
10. Roll the edge of the thumb hole back to clear the thenar eminence.
11. Center the splint along the forearm. Fold over the distal edge to clear the distal palmar crease and continue molding the splint.

### Attaching the metal housings

12. Place one metal housing on the wrist splint at the level of the distal radius and ulna and mark where the housing will attach (see photo A).
13. If necessary, bend the housing slightly so it lies flush against the splint as follows:
  - To prevent damaging the three screws that are protruding from the housing, use the Allen key to screw them into the housing.
  - Use two crescent wrenches and manually bend the housing, or use a vise as follows:
    - Open the jaws 1 1/2" (3.5 cm).
    - Place the housing between the jaws so the ends rest on top of the edges of the jaws and the tubing hole fits between the jaws.
    - From above, place a 3/8" (1.9cm) dowel horizontally between the screw holes and hit it with a hammer to bend the housing.

14. Position the second housing on the lateral border of the forearm section of the elbow splint just distal to the elbow crease and mark where the housing will attach (B).
15. Punch or drill 1/4" (6.4mm) holes in both splints.
16. Attach the housings with the 2-part screws by putting the outer, female portion of the screw through the splint and the male portion through the housing.

To attach the tubing and complete the splint, see instructions 20 through 35.

## ROLYAN PREFORMED DYNAMIC PRONATION/SUPINATION SPLINT

### Completing the elbow and wrist components

17. Check the elbow and wrist components for proper fit. The two splints should be at least 1 1/2" (3.8 cm) apart on the forearm. To trim or modify the splints, spot heat them with a heat gun or in a heat pan of water at a temperature of 150° to 160°F (65°–70°C) until pliable. **Heating the entire splint may cause it to lose its preformed shape. Polyform has a high degree of drape—for better results when modifying the material, shape it with stroking movements.**
18. If desired, apply pre-cut padding to the wrist component.
19. Attach the metal housing to the elbow component using the pre-drilled holes that will be on the lateral border of the forearm. Attach the housing with the 2-part screws by putting the outer, female portion of the screw through the splint and the male portion through the housing.

## **BOTH PREFORMED SPLINT AND KIT**

### **Attaching the tubing**

20. Loosen the three screws in the housings to allow the tubing to fit into the housings. Place one end of the tubing into the housing on the wrist component so the end of the tubing lies flush with the distal end of the housing. Tighten the screws on the housing securely around the tubing (C).
21. Position the other end of the tubing in the housing on the elbow splint, but leave it loose for future adjustments (C).

22. Apply straps as appropriate. On the elbow splint, the distal strap should be distal to the elbow crease but proximal to the housing.
- Preformed Splint:** The straight strap and hook piece included with the wrist component should be used to secure the distal portion of the wrist component.

### **Applying the splint to increase forearm supination**

23. With the elbow splint on the patient, align the wrist splint with the hand so it is parallel to the extremity (D).

24. Bring the wrist splint and tubing over the radial border of the forearm so the tubing crosses over the dorsal surface of the arm (E). Apply the splint.
25. Adjust the length of the tubing by pulling it through the elbow housing. Allow enough space between the tubing and the forearm to clear one finger's width.
26. Tighten the elbow splint housing around the tubing.
27. To attain the appropriate level of force, remove the wrist splint and align it parallel to the dorsum of the forearm (F). Rotate the splint and tubing initially 2½ turns in pronation.
28. See instructions 34 and 35 for final splint adjustments.

### **Applying the splint to increase forearm pronation**

29. With the elbow splint on the patient, align the wrist splint with the hand so it is parallel to the extremity (D).
30. Bring the wrist splint and tubing under the ulnar border of the forearm so the tubing lies along the ulnar border of the arm (G). Apply the splint.

31. Adjust the length of the tubing by pulling it through the elbow housing. Allow enough space between the tubing and the forearm to clear one finger's width.
32. Tighten the elbow splint housing around the tubing.
33. To attain the appropriate level of force, remove the wrist splint and align it parallel to the dorsum of the forearm (H). Rotate the splint and tubing 2 turns in supination.

### **Making final adjustments**

34. Apply the splint. The force of traction can be increased or decreased by increasing or decreasing the number of rotations or by loosening the screws in the proximal metal housing and shortening the tubing.
35. Excess tubing can be trimmed, but allow at least 1½" (3.8 cm) of tubing to extend past the housing. If the splint will be used to increase both pronation and supination, trim the excess tubing when the splint is positioned for supination because of the extra length needed.

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**Patterson Medical, Ltd.**