PGA

Instructions for Use

POLYGLYCOLIC ACID (ABSORBABLE SURGICAL SUTURES, USP)

DESCRIPTION

Polyglycolic Acid is a synthetic absorbable sterile surgical suture composed of homopolymers, of glycolide or copolymers made from 90% glycolide and 10% L-lactide.

The braided suture is dyed violet, green or undyed natural beige. Polyglycolic Acid meets requirements established by the United States Pharmacopeia (USP) for Synthetic Absorbable Surgical Suture.

INDICATIONS

Polyglycolic Acid is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures, but not for use in cardiovascular and neurological procedures.

ACTIONS

Polyglycolic Acid elicits a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of Polyglycolic Acid synthetic absorbable sutures occurs by means of hydrolysis, where the polymer degrades to either glycolic acid or glycolic and lactic acids which are subsequently absorbed and metabolized by the body. Absorption begins and as a loss of tensile strength with appreciable loss of mass. Implantation studies in animals indicate that Polyglycolic Acid retains approximately 65% of its original tensile strength at two weeks post implantation, with approximately 30% remaining at three weeks. Absorption of Polyglycolic Acid absorbable synthetic suture is essentially complete between 77 and 185 days.

CONTRAINDICATIONS

This suture, being absorbable, should not be used where extended approximation of tissue is required.

WARNINGS

Prolonged contact of this or any other suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Under some circumstances, notably orthopedic procedures, immobilization by external support may be employed at the discretion of the surgeon. Do not resterilize. Discard open unused sutures. Store at room temperature. Avoid prolonged exposure to elevated temperatures.

Users should be familiar with surgical procedures and techniques involving absorbable suture before employing Polyglycolic Acid synthetic absorbable suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.

Since any foreign material in the presence of bacterial contamination may enhance bacterial infectivity, acceptable surgical practice should be followed with respect to drainage and closure of contaminated or infected wounds. The use of this suture may be inappropriate in elderly, malnourished, or debilitated patients, or in patients suffering from cancer, anemia, obesity, diabetes, infection, or other conditions which may delay wound healing. As this is an absorbable suture material, the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of the abdomen, chest, joints, or other sites which may undergo expansion, stretching, or distention, or which may require additional support.

PRECAUTIONS

Skin sutures which must remain in place longer than 7 days may cause localized irritation and should be snipped off or removed as indicated. In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

ADVERSE REACTIONS

Adverse effects associated with the use of this device include wound dehiscence, failure to provide adequate wound support in closure of the abdomen, chest, joints, and other sites where expansion, stretching or distension occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients or in patients suffering from cancer, anemia, obesity, diabetes, infection or other conditions which may delay wound healing, enhanced bacterial infectivity, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, and pain, edema and erythema at the wound site.

HOW SUPPLIED

Polyglycolic Acid synthetic absorbable suture is available in sizes 8-0 through 2. The braided suture is dyed violet, green or undyed natural beige. The suture is supplied sterile in pre-cut lengths or on ligating reels, non-needled or attached to various needle types, in one-dozen boxes. Caution: National or local regulations may restrict this device to sale, distribution and use by or on the order of a physician, veterinarian or other licensed practitioner.

ProAdvantage by NDC
Made in BUSAN, KOREA

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