

ProSpray C-60™

Concentrated Surface Disinfectant/Cleaner

Description:

ProSpray C-60 is a concentrated surface disinfectant and cleaner in one product, approved for use in medical and dental settings.

Water Based Formula with Non-Alcohol Actives:

- Water based formula preferred for a wide variety of surfaces including vinyl, metals, plastics, tubing and faux leather upholstery.
- Will not stain or bleach hard surfaces.
- Stays wet on surfaces longer for full contact time.
- Use around electronics is consistent with manufacturers' recommendations for non-alcohol, non-bleach formula.

Hard on Germs - Not on Bare Hands!

- So safe it is now EPA approved for use after dilution without gloves (in absence of bio-hazards).

Convenient and Pleasant to Use:

- Deodorizer with light lemon scent.
- Reduce inventory: One product cleans and disinfects.
- Formulated for effective action in hard water.
- Approved for 60 day use life after dilution.

Residual Activity:

- Sustained antimicrobial action after application which is not found with alcohol or bleach compounds.

Economical Choice:

- 32 ounces makes 8 gallons of diluted solution.
- Easy consistent dilution: 1/2 ounce unit dose tube makes 16 ounces of diluted solution.
- Concentrate reduces shipping costs and storage space.



Concentrated format saves money and storage space.

Product Specifications

Active Agents	9% o-phenylphenol 1% o-benzyl-p-chlorophenol
Dilution	1:32
EPA Reg. #	46851-1
Solution Color	Natural color
Odor	Alcohol, Phenol and Lemon
Use Life	60 days for diluted unused solution
Shelf Life	2 years from date of manufacture



Description	Case Qty	Order#
32 oz. Tip & Pour Measuring Bottle	12	PSC60/32-1
Operator Pack - (48) 1/2 oz. unit dose tubes	3	PSC60/48-1
Intro Kit - (24) 1/2 oz. unit dose tubes, (2) 16 oz. empty spray bottles, (2) squirt caps	1	PSC60/INT-1
16 oz. Empty Spray Bottle	6	PSC60SB-1
Optional Squirt Cap	6	CAPSQUIRT-1

Broad Spectrum Disinfecting Power in One Minute:

- Viruses including HIV-1 and Herpes simplex virus type 1 & 2.
- Vegetative bacteria including Staphylococcus aureus, Methicillin-resistant Staphylococcus aureus, Salmonella enterica and Pseudomonas aeruginosa.
- Fungi including Trichophyton mentagrophytes.

Inactivates Benchmark M. tuberculosis var. bovis (TB) in Five Minutes:

- Meets CDC guidelines for surface disinfection in the presence of blood and other soils.

Four Applications in One Product:

- Surface cleaner: Contains detergents and surfactants for essential pre-cleaning.
- Intermediate-level disinfectant.
- Ultrasonic instrument cleaning and holding solution.
- Disinfect dental impressions and appliances.

Intended Use:

ProSpray C-60 is a broad spectrum antimicrobial disinfectant and cleaner. It is tuberculocidal and meets CDC guidelines for intermediate-level disinfectant use when blood and body soils are likely to be present. Detergents in the water based formula provide effective cleaning. It may be used on surfaces, as a cleaning solution for instruments and as a disinfectant for dental impressions and appliances.

Product Capabilities:

Kills a broad range of microorganisms including HIV-1, Staphylococcus aureus, Staphylococcus aureus - MRSA HA, Salmonella enterica, Pseudomonas aeruginosa, Herpes simplex virus types 1 & 2 and Trichophyton mentagrophytes (Fungus) in one minute. Kills tuberculosis var. bovis (TB) in five minutes.

Note 1: "Potency against M. tuberculosis has been recognized as a substantial benchmark. Any germicide with a tuberculocidal claim on the label is considered capable of inactivating a broad spectrum of pathogens, including such less resistant organisms as bloodborne pathogens (e.g., Hepatitis B virus, Hepatitis C virus and HIV)."

CDC Guidelines for Infection Control in Dental Healthcare Settings - 2003.
MMWR 2003;52(No.RR-17): p. 63.

Cleans and Disinfects:

Because effective disinfection relies on effective removal of soil, ProSpray C-60 is ideal for both steps - no need for two separate products. Special detergents in ProSpray C-60 clean even dried on blood and other soils.

Safe for Surfaces:

ProSpray C-60 is a water based disinfectant with non-alcohol actives. Water based disinfectant formulas are preferred by upholstery manufacturers. And, unlike alcohol based formulas, ProSpray C-60 will not soften or haze vinyls and plastic. ProSpray C-60 will not corrode metals or stainless steel counters as may occur with bleach and hydrogen peroxide based products. It will not stain or bleach surfaces even with repeated use.

How to Use:

For surface disinfection as a "one-step" cleaner/disinfectant for surfaces that are not visibly soiled:

- Spray or squirt diluted ProSpray C-60 onto surfaces or apply using paper towels or 4 X 4 gauze. Allow surfaces to stay wet for five minutes.

For surface disinfection on surfaces with visible blood and soils:

ProSpray C-60 is an intermediate-level hospital disinfectant with TB kill, meeting CDC Guidelines and OSHA for cleaning and disinfecting surfaces with blood and body soils. Use the 2 step method:

- Clean surfaces by spraying or squirting diluted ProSpray C-60 solution onto surfaces or paper towels or 4 X 4 gauze. Wipe to remove soils. Discard used wipe material. Reapply diluted ProSpray C-60 solution to surfaces and keep wet for a total of five minutes.

For instrument holding and cleaning solution:

- Check instrument manufacturer's instructions. Fill pan or ultrasonic tank with sufficient diluted solution to cover instruments. Soak for 5 minutes. Avoid prolonged soaking over one hour.

For disinfecting dental impressions and appliances:

- Check dental material manufacturer's instructions. ProSpray C-60 is a water based disinfectant recommended to disinfect many types of dental impression materials. It is bleach-free and will not corrode the metal portions of appliances. Rinse and remove gross debris from item. Place impression or appliance inside baggie, thoroughly wet all surfaces with diluted ProSpray C-60 solution, allow 5 minute contact time. Remove and rinse item gently.

MicroCare Medical
6120 East 58th Avenue
Commerce City, CO 80022
(800) 843-3343



Visit Us:
[MicroCare.com/Certol](https://www.microcare.com/Certol)

How to prepare solution for 16 oz. bottle:



Dispense 1/2 oz. (15 ml.) or 1 unit dose packet of ProSpray C-60 to empty bottle. (60 day shelf life begins)



Add tap water slowly to fill bottle to bottom of cap threads. Screw sprayer head or squirt top on bottle. Gently swirl bottle contents to evenly mix solution prior to use.



Clean and disinfect surfaces using spray nozzle or optional squirt cap. Make sure the adjustable spray nozzle is in open position prior to use.

Use Precautions:

A slight soap film may occur upon repeated application of ProSpray C-60 due to detergents in the formula. Simply wipe surfaces off at least weekly with plain water and paper towels. Rinse garments or surfaces where ProSpray C-60 is used prior to application of any bleach product. Store solution away from direct heat and direct sunlight. Avoid freezing during shipping or storage.

Availability:

Stocked by MicroCare Medical distributors. Contact MicroCare Medical for a distributor near you or visit MicroCare.com/Certol.

Also Available from MicroCare Medical:



ProSpray™
Ready-to-Use Surface Disinfectant/Cleaner



ProSpray™ wipes
Surface Disinfectant/Cleaner Towelettes

MicroCare Medical
6120 East 58th Avenue
Commerce City, CO 80022
(800) 843-3343



Visit Us:
MicroCare.com/Certol