

School Nurse Ready Room



CLINTON
Ready Room™

an all-inclusive Clinton exclusive

The best and most affordable way to quickly furnish and entire school nurse's office. What you see is what you get: our most popular recovery couch; 48" wall and base cabinets with stainless steel sink and faucet; single, stainless steel glove box holder; 5-leg pneumatic stool; stainless steel waste receptacle and coordinating padded side chair.
(no model substitutions, your choice of Classic Clinton laminate and upholstery colors)

MODEL 3600-27-RR

PACKAGE INCLUDES:

- 3600-27 Chrome Leg Couch in your choice of upholstery colors
- 8248 Wall Cabinets in your choice of laminate colors
- 8048-022 Base Cabinets in your choice of laminate colors with stainless steel sink and chrome faucet set included (choose right or left location or no sink. If no location is selected, sink will be placed as shown.)
- 21335 in your choice of upholstery colors
- GS-3000 Single, Stainless Steel,* Glove Box Holder
- TR-20S 20-Quart, Medium Round, Stainless Steel,* Waste Receptacle
- C-40B Premium, Padded, Side Chair in your choice of upholstery colors

* Not MRI Compatible

3600-27 U-Leg Couch

- Dual frame construction with separate base & headrest sections
- Non-adjustable wedge headrest
- Black, powder-coated, steel U-legs
- Paper dispenser included
- 2" firm foam padding
- Heavyweight, knitbacked upholstery
- 250 lbs. load capacity under normal use



SPECIFICATIONS

72" L x 27" W x 18" H

Prop 65 Warning



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black

8248 Wall Cabinet

- Wall cabinet with 2 doors and 1 adjustable shelf
- Easy-clean, all laminate surface
- Tough, chip resistant edges
- Unit is pre-assembled to save installation time
- Cabinet easily mounts with hook and rail mounting system
- Satin finish pulls
- Choice of solid color or woodgrain laminates
- Fully adjustable, concealed, soft-close hinges

CLASSIC CLINTON
• LAMINATE CASEWORK •




8048-022 Base Cabinet

- Base cabinet with 2 drawers and 4 doors with 2 adjustable shelves, plus stainless steel sink and chrome faucet set. *(Choose right or left location or no sink. If no location is selected, sink will be placed as shown.)*
- Easy-clean, all laminate surfaces
- Tough, chip resistant edges
- Unit is pre-assembled to save installation time
- Fully adjustable, soft-close hinges
- 3.25" backsplash included with top
- Satin chrome pulls
- Choice of solid color or woodgrain laminates
- Customizable with options



Prop 65 Warning

1  **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

8248 SPECIFICATIONS

48" W x 12" D x 24" H

8048 SPECIFICATIONS

48" W x 18" D x 35" H

Finishes



Maple



Dark Cherry



Gray




Slate Gray

2135 Standard Series, 5-Leg, Pneumatic Stool

- Thick padded seat for comfort
- 24" Black nylon, 5-leg base
- Nylon dual wheel casters
- Easy, pneumatic height adjustment
- Premium, stain-resistant, woven, knit-backed, vinyl upholstery
- Maximum load capacity 250 lbs. under normal use



Prop 65 Warning

2  **WARNING:** This product can expose you to chemicals, including diisononyl phthalate (DINP), which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

SPECIFICATIONS

Seat Diameter 16"

Base Diameter 24"

Height Range 18"-23"

Upholstery Colors



Slate Blue



Wedgewood



Royal Blue



Desert Tan



Burgundy



Country Mist



Warm Gray



Gunmetal



Black

GS-3000 Single, Stainless Steel,* Glove Box Holder

- Stainless steel construction
- Open design allows for easy glove identification
- All units mount horizontal or vertical
- Mounting hardware included
- Universal size fits most glove boxes

*Clinton stainless steel items are not M.R.I. compatible.



Shown in the horizontal position

SPECIFICATIONS

10" L x 3.75" D x 5.5" H

Prop 65 Warning

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

TR-20S Medium, Round, Stainless Steel* Waste Receptacle

- All stainless steel outer shell and lid
- Step-on foot pedal operated lid
- One-piece, plastic base keeps metal from contact with floor surface
- Chrome Handle make moving and cleaning easy
- 2 warning labels included
- 20 L capacity (21.13 quarts)

*Stainless steel waste receptacles are not MRI compatible



Prop 65 Warning

WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1

SPECIFICATIONS

10" D X 24"H

C-40B Premium, Padded, Black Frame Side Chair Without Arms

- Black powder coated frame finish
- Heavy duty 1-1/2" diameter legs
- Solid plywood seat and back
- Wall saver guard
- 2" thick, firm foam padded seat and back
- Contoured seat and back
- Adjustable leg levers
- Premium, knit-backed upholstery
- 300 lbs. load capacity under normal use



Prop 65 Warning

WARNING: This product can expose you to chemicals including 1, 3-Butadiene which is known to the State of California to cause cancer and reproductive harm. For more information go to: www.P65Warning.ca.gov.

7

SPECIFICATIONS

Seat 18" W x 18" D

Seat Height 18"

Upholstery Colors



Slate Blue

Wedgewood

Royal Blue

Desert Tan

Burgundy

Country Mist

Warm Gray

Gunmetal

Black