

Re-Order Number	Description	Size
757108	Gel-Foam - Soft Base	16" x 16"
757110	Foam Only - Soft Base	18" x 16"
757112	Foam only, Soft Base	20" x 16"
757115	With Gel Pad, Soft Base	18" x 16"
757117	Gel-Foam - Soft Base	20" x 16"
757120	Foam only, Firm Base	18" x 16"
757122	Foam only, Firm Base	20" x 16"
757125	With Gel Pad, Firm Base	18" x 16"
757127	Gel-Foam, Firm Base	20" x 16"
757130	Gel-Foam Cushion for Geri-Chair	20" x 18"



### Purpose

Provides restraint-free positioning for wheelchair users with forward thrust and sacral sitting problems. Maintains pelvic neutrality. Designed as a comfortable alternative to pommel-style cushions, especially for males.

### Important Notice

This cushion has tie straps and should be used at all times. Although this cushion is designed to help prevent the resident from sliding off a chair, the resident should be monitored as needed. If the resident should be monitored as

### Easy Care

The outer cover can be wiped clean and is removable. It can be hand or machine washed, air or machine dried at temperatures less than 140°F (60°C). The inner cover is heat sealed and can be wiped clean with a standard disinfectants.

