

# Imagination Series, Space Place, Pediatric Scale Table

## **MODEL 7835**


## **TABLE FEATURES**

- Clinton original illustrations
- Easy to clean laminate surfaces
- All graphics are permanently embedded in the high pressure laminate and will not peel off
- Graphics are on one long side and two short ends. Back is slate gray. Only manufactured as pictured.
- Comfortable work height of 36"
- Top is rimmed on three sides for safety
- Top's front and back rim edges are beveled and the corners rounded for safety
- Easy-clean, replaceable center pad
- Built-in infantometer with folding ends for measuring height
- Easy access to weight, measurement and storage
- Choose to access concealed paper dispenser from either end
- Recessed front toe kick
- Easy-clean all laminate base
- Ample 2 door and 1 drawer storage
- Soft-close, door hinges

MODEL 7835 SPECIFICATIONS

43"L x 24"W x 38.5"H

# California residents please be advised of the following, pursuant to Proposition 65:







- Dependable, accurate Doran® scale
- Easy-read digital display
- Auto calibration on start-up
- Auto shut-off
- Electric or battery operation (6 AA batteries required)
- Weights in increments of lbs., ozs. or kilos
- Minimum registered weight is 5 lbs.—Maximum weight for scale is 40 lbs.

Options



066 Door & Inside Latch Combo



055 Drawer Locks

Upholstery



## Features



Infantometer Built-in infantometer with folding ends for measuring height



Paper Dispenser Choose to access concealed paper dispenser from either end



**Rimmed Top with Beveled Edges** Top is rimmed on three sides for safety and the front and back rim edges are beveled and the corners rounded for safety



Padded Top Easy-clean, padded, upholstered top for comfort



**Scale Head** Easy-read, digital display