

# Secure II®

## Med/Surg Bed

**Available with Stryker's patented Chaperone® bed exit system and Zone Control® technology**

**Single-button, in-bed scale provides accurate, repeatable weights with the patient in any position**

**Open architecture allows the flexibility to choose any sleep surface**

**Retractable bed frame helps keep patients close to their belongings regardless of bed articulation**

**Fall prevention educational program and CEU course available**

The combination of Stryker's Chaperone bed exit system and Fall Prevention Program, coupled with the feature of bed retractability, form the most comprehensive fall prevention solution available.



### The Standard in Patient Care

Responsive engineering over time has resulted in a med/surg bed that offers the highest level of security, intuitive operation and flexible functionality. Fully-customizable through a host of features designed around patient and caregiver satisfaction, the Secure II represents the benchmark in patient handling.



Optional Zoom® drive system minimizes the effort needed for patient transport while encouraging proper ergonomic posture.



Intermediate siderail position and low bed height assist with patient entry and exit.



Chaperone accurately senses body positioning by constantly tracking the patient's center of gravity.

# Secure II<sup>®</sup>

## Med/Surg Bed

**stryker<sup>®</sup>**



### Medical

6300 South Sprinkle Road  
Kalamazoo, Mich. 49001 USA  
t: 269 329 2100 f: 269 329 2311  
toll free: 800 787 9537

www.stryker.com

### Standard Features

- 6" (15 cm) casters
- Dual pedestal design
- CPR release
- In-rail fowler/knee gatch control
- Electronic function lock-out control
- Four-wheel, steel-ring brake system with centrally-located activation pedals
- Backrest angle indicator
- Locking-caster steering
- Four traction sockets
- Removable headboard with CPR capability
- Retractable litter top
- Four foley bag hooks
- Four IV receptacles
- Patient restraint locations
- Trend./reverse Trend. controls

### EX Model: All Standard Features Plus

- Chaperone<sup>®</sup> center-of-gravity bed exit system
- Scale with isolated foley bag hooks and Trend. display
- Stryker pendant port

### ZX Model: EX Features Plus

- Zoom<sup>®</sup> Motorized Drive System

### Specifications

• <b>Model Number</b>		
Domestic	3002 (115 VAC)	
International	3221 (230 VAC)	
• <b>Overall Length</b>	93"	(236 cm)
• <b>Overall Width</b>		
Siderails Up	42.5"	(108 cm)
Siderails Down	39"	(99 cm)
• <b>Weight Capacity</b>	500 lbs	(228 kg)
• <b>Height Range (to litter top)</b>		
High	30"	(76 cm)
Low (standard with 6" casters)	16"	(41 cm)
Low (ZX option)	19.75"	(50 cm)
• <b>Litter Positioning</b>		
Backrest	0° - 60°	
Knee Gatch	0° - 40°	
Trend./reverse Trend.	±12° (± 10° with ZX option)	
Retraction	14"	(36 cm)
• <b>Patient Surface</b>	35" x 84"	(89 x 213 cm)
• <b>Caster Diameter</b>		
Standard	6"	(15 cm)
Optional	8"	(20 cm)

### Electronics

• <b>Agency Approvals</b>	C-US UL (3002), CE (3221)
• <b>Volts</b>	115 VAC (Domestic) / 230 VAC (International)
• <b>Ampere Rating</b>	7A (Domestic) / 4A (International) 10A (rating for 115 volt outlet option)
• <b>Current Leakage</b>	<100 microamperes
• <b>Frequency</b>	60 HZ (Domestic) / 50-60 HZ (International)
• <b>Hospital Grade Plug</b>	Yes

*Stryker reserves the right to change specifications without notice.*

### Optional Features

- In-bed scale system
- Chaperone center-of-gravity bed exit system
- Chaperone with Zone Control<sup>®</sup>
- Siderail communications—includes nurse call with speaker, TV, radio, volume, room/read lights
- Smart TV—includes closed caption, channel up/down and mute
- Mattresses: management and prevention
- Adapter bracket for patient helper (not available with Zoom)
- Bed extender
- Removable CPR board
- Defibrillator tray
- Pendant control
- Pillow speaker interface
- Pump holder mount on footboard
- Siderail pads
- Transducer mount—siderail or IV pole
- Traction sleeve adapter kits
- Upright oxygen bottle holder
- Wallsaver<sup>™</sup> quick-release connecting kit
- 8" (20 cm) Omni Surface<sup>™</sup> casters with wheel covers
- IV poles
  - Permanent folding two-stage
  - Removable
- 110 volt power outlet (not available with Zoom)
- Prop rod

Copyright © 2005 Stryker  
LSS-453-06C  
PC/SCG-5M-0805  
Mkt Lit-78-02022005 Rev C