

Let's protect what really matters... our people.



Offering the highest quality and best value in personal protection and wound care products

Disposable Gowns

Comfortable long sleeve: disposable, full-length gowns



Isolation Gowns

Latex Free;
Waist and neck ties;
Economical;
Fluid resistant;
One ply material making it cooler to wear;
Sizing slightly larger than most other popular isolation gowns for easier wear and removal during an emergency situation;
Elastic cuffs;

A8009 – Regular, yellow
A8010 – Extra large, yellow
D8009 – Regular, yellow
D8010 – Extra large, yellow
A8009-B-Cuff – Regular, yellow
A8010-B-Cuff-Blue – Extra large, blue
Packaging: 10/bag, 5 bags/case

A69975 – Premium Protection, yellow
A69979 – Premium Protection, yellow
Packaging : 10/bag, 10bags/case

8023 – Regular, blue
A8015 – Regular, pink
Packaging: 10/bag, 5bags/case



Lab Coats

Latex Free;
Premium;

Thick and durable three-ply material for better protection and comfort;
Snap front buttons;
Knitted cuffs and collar;
2 front pockets;
Highly fluid repellant;

8016 – Small, sky blue
8017 – Medium, sky blue
8018 – Large, sky blue
8019 – Extra large, sky blue

Packaging: 10/bag, 5 bags/case

A8041 – Small, white, no pocket
A8042 – Medium, white, no pocket
A8043 – Large, white, no pocket
A8044 – Extra large, white, no pocket
A8045 – 2X-Large, white, no pocket
A8046 – 3X-Large, white, no pocket
A8048 – 4X-Large, white, no pocket
A8049 – 7X-Large, white, no pocket
A8050 – 9X-Large, white, no pocket

Packaging: 10/bag, 5 bags/case

Impervious Gowns

Latex Free;
Knitted cuffs;
Waist and neck ties;
Superior;
Highly fluid repellant (laminated coating);
Complete body coverage;

8022 – Regular, white
8024 – Extra large, white
Packaging: 10/bag, 5 bags/case

8025 – Regular, blue, with thumbloops
8026 – Extra-large, blue, with thumbloops
A8027 – Premium, regular, blue, with thumbloops
Packaging: 15/box, 5boxes/case

Impervious Gowns with Thumbloops

Shoe Covers

Soft durable and lightweight shoe covers
Premium quality;
Universal size;
Fluid resistant;
Sewn seam offering greater tear strength;
Expandable to male size 13 while still comfortably fitting around a woman size and providing a snug coverage;
Elastic opening;
Competitive price.



Non Conductive

8008 – Blue
8008-XL – Blue, X-Large
Packaging: 100/box, 3 boxes/case

Non Skid, Non Conductive

Non skid tread prints improve traction and reduce chances of slipping.

8007 – Bleu
8006 – Bleu, X-Large
Packaging: 100/box, 3 boxes/case



AMD RITMED INC. CANADA
1200, 55th AVENUE
LACHINE, QUEBEC
CANADA, H8T 3J8



TEL: 888-441-2120 FAX: 888-434-4988
www.amdritmed.com

Bouffant Caps

Cool and functional covers designed for hygiene and protection, where complete coverage is required.



Latex free,
Premium quality,
Fluid resistant,
Soft elastic headband for maximum coverage,
Spunbond polypropylene construction for
durability, ventilation and comfort.

Bouffant Caps

8000 – 21", white
8001 – 21", blue
8003 – 24", white
8004 – 24", blue

Packaging: 100/box, 5 boxes/case

Masks

The perfect combination of safety, quality and comfort!

Our mask products exceed industry standards for filtration and breathability. They contain full flexible nosepieces with shingle pleats, have latex free earloops and fluid resistant inner and outer material.

DISTECH

Latex free masks with earloop

2115 – Blue
2115K – White, small
2116 – Pink
2117 – Yellow

Packaging: 50/box, 10 boxes/case



N95

Certified N95 NIOSH Respirator approved
masks

A2320 – Duck Bill
M2321 – Cone

Packaging: 20/box, 8 boxes/case
Latex free and good breathability



PREMIER

Ultra fluid resistant inner and outer layers.

2015 - Blue
2016 – Pink
2017 – Yellow
2018 – Teal

Also available in Fog Free

Packaging: 50/box, 10 boxes/case




CONE FACE MASK

Specially designed smooth edges for a soft comfortable
fit. Featuring a comfortable and practical elastic
headband.

2020 – Blue
2021 – Pink

Packaging: 50/box, 12 boxes/case



PRO-SHIELD

Fluid resistant. Full coverage pre-attached
wraparound and optically clear fog free visor.

2125 – Teal, earloop
2126 – Blue, tie-on

Packaging: 25/box, 4 boxes/case



TIE-ON

Fluid resistant outer layer and a specially treated
tissue inner layer. Sonically sealed ties on outside

2000 – Blue
Also available in Fog Free

Packaging: 50/box, 6 boxes/case



2900 – Protective eyewear

Packaging: 10/box, 6 boxes/case

