

## Mobile Aneroid Sphygmomanometer

- Height adjustable aneroid sphygmomanometer
- Luminescent dial
- Large, 6", easy-to-read gauge mounted on a mobile stand
- Gauge mounted on a 110° swivel range bracket
- Accumax™ nylon cuff with 8' coiled tubing
- Storage basket attached on stand
- One (1) year warranty for manufacturing defects only

**V223** Adult, Latex Free Bladder Black

**223B** Adult, Latex Bladder Black

**Caution: This product contains natural rubber latex, which may cause allergic reactions.**

Cuff Class	Limb Range (Inches)	Cuff Dimensions (Inches)*
Newborn	3.75" – 5.25"	7" W x 2" H
Infant	5" – 7.5"	10.5" W x 3" H
Child	7" – 10.5"	13.5" W x 4" H
Adult	9.5" – 15"	20.5" W x 5.5" H
LG Adult	13" – 19"	25" W x 6.5" H
Thigh	17" – 26"	29.5" W x 8.5" H

\* These are subjective measurements and can vary  $\pm 0.25"$

