

Adjustable Forearm Platforms

HCPCS: E1399



- Height and depth adjustable
- Rotates inward and outward
- Handgrips extend forward or backward, supinate, pronate and are angle adjustable
- Contoured armrest is molded into one piece with a flexible edge for added comfort
- Hook and loop fastener straps prevent arms from slipping
- May be mounted anywhere along the handlebar, including the front Position is determined by patients weight bearing needs
- May be mounted facing toward the inside or outside of the handlebar for lateral adjustability
- Each side is individually adjustable

PRODUCTS	Description
CE 1035 FP	rtl