

Bariatric Offset-Handle Canes

HCPCS: E0100



FIGURE A




FIGURE B

- Manufactured with steel tubing
- Easy-to-use, dual-button height adjustment with locking ring Locking ring prevents rattling (Figure A)
- Comfortable vinyl hand grip with wrist strap
- Model 10305 Weight Capacity 500 lbs
- Model 10318 Weight Capacity 600 lbs

PRODUCTS

Description

10305-6	Cane,HD,Offset,Black,6/cs
10318-6	Offset Cane,HD,Extra Long,6/cs

SPECIFICATIONS

Description

Handle Height	29.75" - 37.75" 37" - 46"
Warranty	Limited Lifetime
Weight Capacity	600 lbs