




SUIT UP

Medicom has you covered

Medicom protective apparel is made of high-quality materials and designed to provide complete coverage that complies with all applicable regulations. It provides the reliable protection you need without compromising comfort or style and helps create a reassuring environment for patients.

Benefits of single-use protective apparel



Reduced risk of cross-contamination



Always clean, as packaging remains sealed until use



No laundering required

Protective Apparel Best Practices

Protective apparel must be worn and disposed of correctly. We encourage you to share these guidelines with your staff:

Donning protective apparel

- Perform hand hygiene before removing item from package.
- Select appropriate garment for intended use (length, coverage, fluid resistance, etc.)
- Fasten correctly: tied at neck and waist to ensure comfortable fit, adequate coverage and freedom of movement.
- Protective apparel should not be worn in restrooms, personal offices or break rooms.

Doffing protective apparel

- Remove apparel immediately after use. Do not re-use.
- Perform proper hand hygiene after removing apparel to prevent cross-contamination.
- Turn garment inside out when removing to keep infectious material inside.
- Designate a safe receptacle in which to dispose of protective apparel after use.

When to change apparel

Change protective apparel between patients and as soon as possible when penetrated by blood or other potentially infectious fluids.

SafeWear™

Premium Protective Apparel



Discover the complete range

Available in Hipster Jackets, High Performance Lab Coats and Form-Fit Isolation Gowns, made from top-quality materials and designed to provide reliable barrier protection without sacrificing comfort or style.

FIT

- Thoughtful designs offer ease of movement, as well as secure closure and easy removal.
- Special features such as front snaps, back slits and neck and waist ties, provide a safe, comfortable and personalized fit.
- Lightweight, soft, supple fabrics keep wearers cool and comfortable.
- SafeWear protective apparel is latex free.
- Available in a wide range of sizes, from XS - 3XL.

FUNCTION

- SafeWear protective apparel complies with recommended OSHA best practices.
- Available in a wide range of material weights and polypropylene/SMS variants to meet a variety of needs while maintaining breathability.
- Premium materials provide superior protection, reliability and durability.
- Hipster Jackets and Lab Coats are made from AAMI PB70 Level 2 tested material.

Hipster Jacket™ Premium hip length jacket

Item No.	Description	Unit/Pkg	Pkg/Case
8101*	Deep Blue	12	5
8102	Pretty Pink	12	5
8103*	Plum Purple	12	5
8104*	Soft Blue	12	5
8105*	White Frost	12	5
8115	Tropical Teal	12	5
8116	Poppy Pink	12	5
8118	Jet Black	12	5

Form-Fit Isolation Gown™ Premium isolation gown

Item No.	Description	Unit/Pkg	Pkg/Case
8113*	Sunny Yellow	12	5
8114*	Bright Blue	12	5

High Performance Lab Coat™ Premium lab coat

Item No.	Description	Unit/Pkg	Pkg/Case
8107*	Deep Blue	12	5
8108	Pretty Pink	12	5
8109*	Plum Purple	12	5
8110*	White Frost	12	5
8112*	Soft Blue	12	5
8117	Tropical Teal	12	5
8119	Jet Black	12	5

*All colours available in the USA. Colours with an asterisk available in Canada.

Size:

XS: Extra Small A: Small B: Medium C: Large D: Extra Large E: 2x Large F: 3X Large

Hipster Jackets and High Performance Lab Coats: All colours available in Small, Medium, Large and Extra Large. 2X Large and 3X Large are available in Deep Blue (3X Large available only in the USA). Extra small is available in Plum Purple for Hipster Jackets only.

SafeBasics™

SafeBasics protective apparel provides a basic level of protection at an attractive price.

Isolation Gowns

Single-use gown for procedures requiring protection from spray and splatter. Fluid-repellent design provides added protection against cross-contamination.



Lab Gowns

Front snap closure, knitted cuffs and collar for comfortable, easy on-and-off fit.



Fluid Impervious Gowns

Polyethylene coating provides fluid resistance for procedures involving the potential for cross-contamination.



Item No.	Description	Unit/Pkg	Pkg/Case
8009	Isolation Gown Reg - Yellow	10	5
8010	Isolation Gown XL - Yellow	10	5
8023	Isolation Gown Reg - Blue	10	5

Item No.	Description	Unit/Pkg	Pkg/Case
8016	Lab Gown S - Sky Blue	10	5
8017	Lab Gown M - Sky Blue	10	5
8018	Lab Gown L - Sky Blue	10	5
8019	Lab Gown XL - Sky Blue	10	5

Item No.	Description	Unit/Pkg	Pkg/Case
8022	Fluid Impervious Gown Reg - White	10	5
8024	Fluid Impervious Gown XL - White	10	5



Bouffant Caps

Hair and head fluid protection with breathable, low-linting design.

Beard Covers

Face protection with an elastic edge and headband loop for optimal fit and comfort.

Item No.	Description	Unit/Pkg	Pkg/Case
8000	Bouffant Caps 21" White	100	5
8001	Bouffant Caps 21" Blue	100	5
8003	Bouffant Caps 24" White	100	5
8004	Bouffant Caps 24" Blue	100	5
8035	Beard Covers, White	100	10

Shoe Covers

Footwear protection with fluid repellent, slip-on covers.



Item No.	Description	Unit/Pkg	Pkg/Case
8006	Shoe Cover Non-Conductive Non-Skid-XL - Blue	100	3
8007	Shoe Cover Non-Conductive Non-Skid-Reg - Blue	100	3
8008	Shoe Cover Non-Conductive Reg - Blue	100	3
8036	Shoe Cover Non-Conductive Non-Skid-XL - White	100	3