

the Landing Strip®

from Comfortex®

Bonded Top Cover

Transfers impact forces outward

High Density Composite Core

Absorbs impact, slows forces

Non-slip Base

Stays where needed

Tapered to the floor Perimeter

Staff and Resident Safety

Sealed Edges

Aid with Infection control



Sizes available
24" or 36"

Nite-Glo™ Option

Nursing is programmed to pick up things on the patient room floor. But the Landing Strip is designed to remain alongside your patient's bed day and night. It guards against injury should a fall occur. With the addition of the Landing Strip Nite-Glo™, the border of the mat stands out, even in a dimly lit room, announcing its presence to both patients and staff.

*It is a thin line between functional and
fantastic. We defined that line...*

Foldable mats



Thick, soft mats




Engineered Excellence




How it works

1. The Landing Strip lies unobtrusively alongside the bed, ready to absorb and redirect impact forces should a fall occur.
2. The combination of firm, compressed foams, bonded top fabric and free-floating bottom fabric react the moment of impact.
3. Deceleration begins and the high density core prevents a sudden bottoming-out against the floor. Like the airbags used by stunt people, downward force can dissipate along the underside of the mat, away from the point of impact preventing a sudden reversal of energy. Softer mats may be able to bounce a falling egg, but can fail under the force generated by a falling person.



A Costly Problem

With rising acuity levels, the risks of patient falls, side-rail entrapment and serious injury are also on the rise. Fall related injuries break bones, spirits and budgets. They are physically and emotionally traumatic to the patient, family members and healthcare providers. The cost of a non-fatal fall averages more than \$20,000 Nationwide. That's more than \$20 billion annually.



Landing Strip™				
24" x 70"		36" x 70"		
Each	Box of 4	Each	Box of 3	Color
0815-91-213	0815-91-338	0815-91-239	0815-91-361	Gray
0815-91-221	0815-91-346	0815-91-247	0815-91-379	Burgundy



www.pattersonmedical.com

Call toll free:
1-800-323-5547
 In Canada call 1-800-665-9200

Fax toll free:
1-800-547-4333
 In Canada call 1-877-566-8611