

Med-Aire Melody Alternating Pressure and Low Air Loss Mattress Replacement System

HCPCS:



FIGURE A

- Active mattress replacement system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and microclimate control
- Designed to prevent, treat and heal pressure ulcers in the home or long term care setting
- Control unit offers 10 minute cycles when alternating and can also be used in static mode
- Adjustable patient weight settings allow for optimal immersion and patient comfort & compliance
- Visual alarm alerts when air pressure is low
- System includes 20 individual air cells to maximize immersion and envelopment with pillow feature for head stability & comfort
- 8" deep air cells are easily removed for convenient servicing
- Removable Nylon cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- CPR valve allows for rapid deflation
- Tubing quickly and easily disconnects with EZ Lock/Quick Release design
- Also available with an audible low pressure alarm (#14026-A)

PRODUCTS

Description

14026	Med Aire Melody Alt Press Sys 35Wx80Lx8H
14026-A	Med Aire Melody w/alarm Sys 35Wx80Lx6H

SPECIFICATIONS

Description

Weight Capacity	350 lbs
-----------------	---------

ACCESSORIES

Description

14427QF	Quick Fill Pump w/adapters, 1 EA/12 CV, EA
---------	--