

SafeBasics® Dry-Back Bibs



TECHNICAL DATA SHEET

FEATURES & BENEFITS

3-Ply Bibs

- Suitable for minor procedures like cleanings, where extra protection is needed.
- Superior absorption and clothing protection with additional layer of paper and plastic backing.
- Excellent protection thanks to horizontal embossed pattern and a unique water-repellent edge.
- Available in a variety of colours to coordinate with Medicom SafeBasics® plastic cups.

RECOMMENDED USE

- Bibs are made of paper intended to absorb saliva or moisture from the oral cavity during procedures. Latex-free.



CERTIFICATION & STANDARDS

ISO 13485 certified manufacturer.
Health Canada registered medical device manufacturer.

STORAGE CONDITIONS

Bibs should be disposed of after each use.
Store in a cool, dry place.



PRODUCT DETAILS

Product Name:	Medicom SafeBasics® Dry-Back Bibs
Product Type:	Single-use, non-sterile, 3 ply tissue + 1 ply poly
Size:	13" x 18" (33 x 45,7 cm)
Quantity per Case:	500 Bibs
Color:	See Appendix 1
Country of Origin:	USA

BIB SPECIFICATIONS

CHARACTERISTIC	SPECIFICATION		
Material Composition	Tissue Paper & Polyethylene Film		
Description	Folded piece of paper with a plastic backing. Embossed horizontal pattern on the body of the bib and vertical embossed patten on the vertical borders.		
Dimension	Unfolded Size	Unfolded Width	Unfolded Length
	13" x 18" (33 x 45,7 cm)	13"	18"
	Folded Size	Folded Width	Folded Length
	4 5/8" x 6 5/8" (11,8 x 16,8 cm)	4 5/8"	6 5/8"

Appendix 1 (List of items)

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment	Marketed in USA as	Marketed in Canada as
8380	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Lavender	3 Ply Tissue + 1 Ply Poly	Non-Sterile	General Purpose	General Purpose
8381	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Silver Grey	3 Ply Tissue + 1 Ply Poly	Non-Sterile	General Purpose	General Purpose
8382	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Blue	3 Ply Tissue + 1 Ply Poly	Non-Sterile	General Purpose	General Purpose
8383	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	White	3 Ply Tissue + 1 Ply Poly	Non-Sterile	General Purpose	General Purpose
8384	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Yellow	3 Ply Tissue + 1 Ply Poly	Non-Sterile	General Purpose	General Purpose
8385	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Green	3 Ply Tissue + 1 Ply Poly	Non-Sterile	General Purpose	General Purpose

Appendix 1 (List of items)

Item Code	Description	Packaging Information	Size	Color	Specific Attribute	Treatment	Marketed in USA as	Marketed in Canada as
8386	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Dusty Rose	3 Ply Tissue + 1 Ply Poly	Non-Sterile	General Purpose	General Purpose
8387	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Beige	3 Ply Tissue + 1 Ply Poly	Non-Sterile	General Purpose	General Purpose
8388	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Peach	3 Ply Tissue + 1 Ply Poly	Non-Sterile	General Purpose	General Purpose
8389	Medicom SafeBasics Dry-Back Bibs	500 Bibs/Case	13" x 18" (33 x 45,7 cm)	Aqua	3 Ply Tissue + 1 Ply Poly	Non-Sterile	General Purpose	General Purpose

AMD Medicom Inc. | 2555 Chemin de l'Aviation, Montreal, Quebec, Canada H9P 2Z2