

R1 shock-proof®

Exact precision that lasts! Made in Germany.



product
design
award

2010



Handheld, portable, aneroid Sphygmomanometer with patent pending shock proof technology. Validated by the British and Irish Hypertension Society. Available in white or black.

- › Patented (pending) shock-proof technology.
- › Special air duct for stable pointer movements.
- › Fast, ergonomic inflation due to optimized bulb design with integrated spoon.
- › Non-resistant opening of air release valve ensures precise and comfortable operation.
- › One piece cuff, washable and disinfectable.
- › 5 years warranty on calibration.
- › Patented (pending) shock-proof technology to height of 120cm / 4 feet
- › Two-component ABS casing.
- › Focus Green Silver 2008 Design Award.
- › Scale: Ø 52 mm white with blue pointer or black with red fluorescent pointer.
- › Valve base and air duct with metal precision valve.
- › Metal threaded tube connection at the top of the casing for fast cuff exchange and ergonomical operation.
- › Maximum error tolerance of +/- 3 mm Hg.
- › Extensive selection of various cuff sizes.
- › Incl. Vinyl zippered bag.
- › Weight: 173 gram.

No. 1250-100	R1 shock-proof® black manometer, 1-tube, w/o cuff	No. 1251-100	R1 shock-proof® white manometer, 1-tube, w/o cuff
No. 1250-107	R1 shock-proof® black, 1-tube, velcro cuff (adult)	No. 1251-107	R1 shock-proof® white, 1-tube, velcro cuff (adult)
No. 1250-150	R1 shock-proof® black, 1-tube, disinfectable one-piece cuff (adult)	No. 1251-150	R1 shock-proof® white, 1-tube, disinfectable one-piece cuff (adult)

* Further items are available, please contact your Riester distributor for further information.

Please visit riester.de for more information