



# AEROSOL DRAINAGE BAGS

## 750 cc Bag • Accepts Standard 22 mm Corrugated Tubing

Dynarex Resp-O<sub>2</sub>™ Aerosol Drainage Bags are designed to prevent accumulated liquid and condensation from reaching patients' airways during vent/aerosol treatment.

- Traps and contains accumulated liquid and condensation
- Self-locking plastic chain secures bag to stable objects
- Not made with natural rubber latex or DEHP



Y-Adaptor

**REORDER NO. 34200**

**PACKAGE QUANTITY**  
50/Case

T-Adaptor

**REORDER NO. 34201**

**PACKAGE QUANTITY**  
50/Case



WE CARE LIKE FAMILY™

1.888 DYNAREX

DYNAREX.COM

10 GLENSHAW STREET | ORANGEBURG, NY 10962

Contact your local sales representative or call to find your nearest distributor.

©2021 Dynarex Corporation



R210428