



WE CARE LIKE FAMILY™



## ECG Snap Electrodes - Disposable

### Snap-style electrodes allow for fast and easy testing

- For ECG/EKG monitoring or diagnosis with the connected equipment as medical sensors
- Ag/AgCl coating on electrode surface for effective signal transmission
- Foam dimensions 50 mm x 55 mm
- ECG Snap Electrodes are for adult use only



Reorder No.	Description	Quantity
7101	50 mm x 55 mm - Adult	200/3/Cs
7102	50 mm x 55 mm - Adult	150/4/Cs
7103	50 mm x 55 mm - Adult	12/50/Cs

See our website for our Table Paper (Reorder No. 4479-4482, 4485-4487)



Find us on:



**dynarex**® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • [www.dynarex.com](http://www.dynarex.com)