

TECHNICAL DATA SHEET
Gauze Sponge
Intended Use

A gauze sponge for external use, device intended for medical purposes, such as to be placed directly on a patient's wound to absorb exudates.

Characteristic		Specification		
Material Composition		100% cotton, Latex-Free		
Physical Properties	Acidity & Alkalinity	Neutral		
	Absorbency	Retention	≥ 9 g/g	
		Submersion	≤ 30 s	
Description		Folded gauze		
Color		White		
Dimension		Size	Length (cm)	Width (cm)
		2" x 2"	5,0	5,0
		3" x 3"	7,5	7,5
		4" x 4"	10,0	10,0
		8" x 4"	19,5	10,0

Additional Information

Storage Condition	Country of Origin
Keep in a cool dry place	China

Appendix 1 List of Item Codes and Attributes

Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Treatment
2120	Medicom SafeBasics All Gauze Cotton Sponges	2 Sponges/Pouch, 50 Pouches/Tray, 30 Trays/Case, 3000 Sponges/Case	2" x 2"	Sterile
2212	Medicom SafeBasics All Gauze Cotton Sponges	1 Sponge/Pouch, 100 Pouches/Tray, 24 Trays/Case, 2400 Sponges/Case	2" x 2"	Sterile
2800	Medicom SafeBasics All Gauze Cotton Sponges	2 Sponges/Pouch, 50 Pouches/Tray, 30 Trays/Case, 3000 Sponges/Case	2" x 2"	Sterile
3001	Medicom SafeBasics All Gauze Cotton Sponges	200 Sponges/Bag, 25 Bags/Case, 5000 Sponges/Case	2" x 2"	Non-Sterile
3002	Medicom SafeBasics All Gauze Cotton Sponges	200 Sponges/Bag, 40 Bags/Case, 8000 Sponges/Case	2" x 2"	Non-Sterile
3003	Medicom SafeBasics All Gauze Cotton Sponges	200 Sponges/Bag, 20 Bags/Case, 4000 Sponges/Case	3" x 3"	Non-Sterile
3004	Medicom SafeBasics All Gauze Cotton Sponges	200 Sponges/Bag, 20 Bags/Case, 4000 Sponges/Case	3" x 3"	Non-Sterile
3005	Medicom SafeBasics All Gauze Cotton Sponges	200 Sponges/Bag, 20 Bags/Case, 4000 Sponges/Case	4" x 4"	Non-Sterile
3006	Medicom SafeBasics All Gauze Cotton Sponges	200 Sponges/Bag, 10 Bags/Case, 2000 Sponges/Case	4" x 4"	Non-Sterile
3009	Medicom SafeBasics All Gauze Cotton Sponges	200 Sponges/Bag, 10 Bags/Case, 2000 Sponges/Case	8" x 4"	Non-Sterile
3120	Medicom SafeBasics All Gauze Cotton Sponges	2 Sponges/Pouch, 40 Pouches/Tray, 30 Trays/Case, 2400 Sponges/Case	3" x 3"	Sterile
3312	Medicom SafeBasics All Gauze Cotton Sponges	1 Sponge/Pouch, 100 Pouches/Tray, 24 Trays/Case, 2400 Sponges/Case	3" x 3"	Sterile
3800	Medicom SafeBasics All Gauze Cotton Sponges	2 Sponges/Pouch, 40 Pouches/Tray, 30 Trays/Case, 2400 Sponges/Case	3" x 3"	Sterile

Item Code	Description	Packaging Information	Size	Treatment
4120	Medicom SafeBasics All Gauze Cotton Sponges	2 Sponges/Pouch, 25 Pouches/Tray, 24 Trays/Case, 1200 Sponges/Case	4" x 4"	Sterile
4412	Medicom SafeBasics All Gauze Cotton Sponges	1 Sponge/Pouch, 100 Pouches/Tray, 12 Trays/Case, 1200 Sponges/Case	4" x 4"	Sterile
4800	Medicom SafeBasics All Gauze Cotton Sponges	2 Sponges/Pouch, 50 Pouches/Tray, 12 Trays/Case, 1200 Sponges/Case	4" x 4"	Sterile