

TECHNICAL DATA SHEET

Preparation Wipes

Intended Use

Spunlace fabric which is soft, highly absorbent and is suitable for a range of hard surface and delicate wiping applications.

Characteristic	Specification		
Material Composition	Rayon and polyester fibers (Latex free).		
Description	White and soft wipes made of non woven fabric manufactured from a blend of cellulosic and synthetic fibers. (Refer to Appendix 1 for specific attribute)		
Color	White		
Ply	4 ply		
Dimension	Size (Open Washcloth)	Width (cm)	Length (cm)
	8" x 13"	20,4	33,0
	10" x 13"	25,4	
	12" x 13"	30,5	

Additional Information

Storage Condition	Country of Origin
Store under ambient temperature and dry conditions	USA

Appendix 1 List of Item Codes and Attributes

Appendix 1 (List of items)

AMD Medicom Inc.

Item Code	Description	Packaging Information	Size	Treatment
SP-40813	AMD-Ritmed Dry Wipes	50 Washcloths/Pack, 10 Packs/Case, 500 Washcloths/Case	8" x 13"	Non-Sterile
SP-41013	AMD-Ritmed Dry Wipes	50 Washcloths/Pack, 10 Packs/Case, 500 Washcloths/Case	10" x 13"	Non-Sterile
SP-41213-1	AMD-Ritmed Dry Wipes	50 Washcloths/Pack, 20 Packs/Case, 1000 Washcloths/Case	12" x 13"	Non-Sterile