

Dedicated to quality. Committed to comfort.

From plunger to tip, UltiCare® safety syringes are precision engineered to provide superior comfort. Each needle is individually inspected using laser technology, ensuring every injection is as comfortable as possible.

UltiCare® safety syringes compared to other brands

	UltiCare	BD Safety-Lok™	Covidien Magellan™
Multi-bevel tip for easier insertion	✓	✓	✓
Micro polished for a smooth finish	✓	✓	✓
Silicone coating for added comfort	✓	✓	✓
Surgical-grade stainless steel needle for superior strength	✓	✓	✓
Bold graduations for accurate medication dosing	✓	✓	✓



UltiCare®

www.ulticare.com

Safety syringes 100-COUNT BOX

low dead space with permanently attached needle

UltiCare safety syringes with permanently attached needles are engineered with features that help reduce the risk of accidental needlestick injuries and minimize medication waste. The permanently attached sheath is preferred by most clinicians and requires minimal change in training.



Insulin Safety Syringes (LOW DEAD SPACE)	NDC#	Product #
1/2cc (50-unit) - 29G x 1/2" (12.7mm)	08222.0325.91	03259
1cc (100-unit) - 29G x 1/2" (12.7mm)	08222.0321.95	03219




Tuberculin Safety Syringes (LOW DEAD SPACE)	NDC#	Product #
1.0mL - 28G x 1/2" (12.7mm)	08222.6300.25	63002
1.0mL - 27G x 5/8" (15.9mm)	08222.6300.01	63000
1.0mL - 27G x 1/2" (12.7mm)	08222.6300.18	63001
1.0mL - 25G x 5/8" (15.9mm)	08222.2515.89	25158
1.0mL - 25G x 1" (25.4mm)	08222.2511.07	25110

Safety syringes 100-COUNT BOX

needle and syringe combo with detachable needle

UltiCare 3mL luer lock safety syringes with detachable needles combine the benefits of luer lock technology with the safety features engineered to reduce the risk of accidental needlestick injuries. The luer lock needles are also color-coded by gauge to help provide quick and easy identification.



3mL Luer Lock Safety Syringes	NDC#	Product #
3.0mL - 25G x 1" (25.4mm)	08222.6300.87	63008
3.0mL - 25G x 5/8" (15.9mm)	08222.6300.70	63007
3.0mL - 23G x 1" (25.4mm)	08222.6300.63	63006
3.0mL - 22G x 1-1/2" (38.1mm)	08222.6300.56	63005
3.0mL - 22G x 1" (25.4mm)	08222.6300.49	63004
3.0mL - 21G x 1-1/2" (38.1mm)	08222.6300.32	63003
3.0mL - Syringe Only	08222.6301.48	63014

