

Med-Aire Plus 8" Alternating Pressure and Low Air Loss Mattress System

HCPCS:



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Active combination therapy system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and microclimate control
- Designed to prevent, treat and heal pressure ulcers in the home or long term care setting
- Digital control unit offers 4 cycle times and adjustable comfort weight settings for optimal immersion and patient satisfaction
- 3 diagnostic audio and visual alarms alert for issues with air pressure, the system, or power
- System includes 20 individual air cells to maximize immersion and envelopment with pillow feature for head stability & comfort
- Advanced Clinical and Care Giver Features: Auto-Firm, Seat Inflate, microprocessing feedback, lock-out and more!
- Cell-on-cell designed air cells prevent "bottoming out" and provide up to 24 hours of power outage protection
- 8" deep air cells are easily removed with quick disconnects for convenient servicing
- Removable 4 way stretch PU cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- CPR valve allows for rapid deflation
- Tubing quickly and easily disconnects with EZ Lock/Quick Release design

PRODUCTS

Description

14029	MedAire8 LowAirLoss.w/Pump system,1/cs
14029-84	Med Aire 8" Low Air Loss APM w/pump 84

SPECIFICATIONS

Description

Control Unit Airflow	8 LPM
Control Unit Alarms	Visual/Audible
Control Unit Cycle Time	10, 15, 20, 25 Minutes
Control Unit Type	Compressor
Mattress Cover Material	Fluid-Resistant Stretch
Warranty	18 Month Limited
Weight Capacity	450 lbs

ACCESSORIES

Description

14427QF	Quick Fill Pump w/adapters, 1 EA/12 CV, EA
14427QF	Quick Fill Pump w/adapters, 1 EA/12 CV, EA

PARTS

Description

14029XP	Digital 3-Outlet Pump 14029/14029DP,1/ea, EA
14029-CELL	Bladder Only,14029,1/ea", EA
14029B	Carry bag for 14029,1/ea, EA

PARTS	Description
14029MN	Mattress Only,14029,3Tube System,1/ea, EA
14013	CPRPlugForLowAirLossMatt,3Port,1ea, EA
14029-CELL-1	Cell Only wLazer Holes,14029,1/ea, EA
DS16C01207	Mattress Hose 3 port, 1/ea, EA
14029-FFTCPLUG	Plug For Connector,14029 Cell,Female,1e, EA
14029C	Cover for 14029,1/ea, EA
15005PC	Power Cord,1/ea, EA
14001-l	L-Shape Connector, LAL,1/ea, EA
14001-t	T-Connector for LAL, 1/ea, EA
14029-FFTC	Connector for 14029 Cell, Female, 1/ea, EA
14028v	OneWayValveForLowAirLossMatt,1/ea, EA
14029-fmtc	Connector for 14029 Cell, Male, 1/ea, EA
14028T	Tubing,7", LAL, 1/ea, EA
14030xp-s	Pump Socket 3 Hose Connection,1/ea, EA
14007-3port	Tubing Adapter 3-Port,14029DP,1/ea, EA
14030hk	Hanging Hooks,14030XP,1pair, PR
14029FN	Filter14029XP,1ea, EA