



Hydraulic Deluxe Chrome-Plated Patient Lift

HCPCS:

- Chrome-plated steel provides maximum strength
- Safely raises or lowers individuals from any stationary position
 6 point Swivel Bar that can accommodate, 2 or 4 sling straps, plus 2 chain connections
- High-performance hydraulics raise or lower individuals gradually and safely
- Easy-to-operate caster brakes provide security
- Includes sling chains

PRODUCTS	Description
13023CSET	CasterSet,13023, 4/bx
SPECIFICATIONS	Description
Warranty	Limited Lifetime (frame), 1 Year Limited (pump)
Weight Capacity	450 lbs
Boom Height	17.5" - 62"
Carton Shipping Weight	81 lbs
Casters	5"
Weight	72 lbs