

## ADA Lap Accessories



ADA SB

### The Lap Shopper

The Lap Shopper provides hands free convenience.

This enables the shopper with limited dexterity to have both hands free to pick up merchandise from the shelf.



### Fishing Rod Holders

Go fresh or salt water fishing with this great rod holder! The rod holder plugs into the extension tube of the ADA Lap Base.

Available in either fresh water or salt water models. For use with an Extension Pole.

Compatible with Scotty Hardware.

ADA FW Fresh Water Rod Holder  
ADA SW Salt Water Rod Holder



### The ADA Lap Food Tray in use by Gerry Price

Before his workplace accident in 1995, Gerry worked in Engineering and Operations. Since becoming a Quadriplegic, he understands that to be alive and healthy the disabled need to get out and be active in their communities.

The ADA Lap™ does this by allowing the wheelchair user to be active and provides them with a sturdy lap mount. This leaves the hands free. The ADA Lap can serve as a:

- Camera Mount
- Fishing Rod Holder
- Umbrella Holder
- Shopping Basket Mount
- Snap Lap Desk
- Food Tray Holder and more!

### More Fine Products from Symmetric Designs

Savant  
Wheelchair  
Headrest



Headmaster  
Collar



Wheelchair  
Headrest  
Hardware



[www.symmetric-designs.com](http://www.symmetric-designs.com) 

# ADA Lap™ Solutions for an Active Lifestyle



## Universal Design to get you going!

Add a lap desk, shopping basket, fishing rod, umbrella, still or video camera, food tray or whatever you choose.

The ADA Lap™ simply slides under your power or manual wheelchair cushion.



Symmetric Designs

125 Knott Place, Salt Spring Island, BC V8K 2M4, CANADA  
email: [sales@symmetric-designs.com](mailto:sales@symmetric-designs.com)  
[www.symmetric-designs.com](http://www.symmetric-designs.com)  
Tel: 800-537-1724 Fax: 250-537-1998

# ADA Lap System

# ADA Lap Accessories

## Universal Design To Get You Going!

The ADA Lap™ simply slides under your power or manual wheelchair cushion to provide a versatile base for attachments.




**ADA Lap** (Pronounced add-a-lap) Solutions provide a universal support for wheelchair users.

The easy-to-use base simply slides under the wheelchair cushion. It is held in place by the weight of the user. Height adjustment is by a quick release clamp and the vertical pole is easily removeable for transfer.



ADA 01 - Base Unit



ADA TR



ADA FT



ADA LD

## Lexan Tray

This strong polycarbonate tray provides a firm transparent surface for writing


## Food Tray

With a lightly textured bottom and curved lip, this 11" x 14" sturdy tray will keep drinks and plates from sliding off.

Also great for containing any small objects that may roll or slide on a flat surface.

## SNAP Lap Desk

A small attache case that opens to contain your papers, pens, etc., yet doubles as a desk. Use with a non-slip surface for keeping control of slippery objects. No more juggling items on your lap!



ADA PL-12  
ADA PL-18



ADA AA



ADA CM

## Extension Poles

Available in either 12" or 18" lengths, extension poles have an insert for camera mounts/tripods, Umbrella holders, Angle Adapters and fishing rod holders.

## Angle Adapter

This versatile adapter comes with a removable top plate surfaced with velcro. It can be used to angle trays and other attachments. For use with an extension pole.

## Camera Mount

The ADA Lap camera mount/tripod allows small height adjustments using a hand crank as well as allowing camera tilt function. This camera mount can be used on the ADA Lap or freestanding.

For use with an extension pole when mounted.

## Umbrella Connector

The ADA Lap Umbrella connector can be screwed into a pre-drilled hole in any wooden handled umbrella. Also available preinstalled with umbrella. For use with an extension pole.



ADA UC or UB