

## 3pp® ThumSling®



The 3pp® ThumSling® is a soft, lightweight splint that applies light compression to relieve arthritis pain at the base or CMC joint of the thumb. The ThumSling gently holds the joint without limiting motion or function. Also used for ligament sprains and thumb pain.

## **Features**

- The wrap-around sling adjusts to the wearer's level of comfortable support.
- Designed to hold the CMC joint in the proper position so the wearer can perform daily activities with less pain or fear of dropping things due to CMC arthritis.
- Easily increase the level of support by adding the flexible stay (included) for additional control wherever it's needed.
- Hand or machine washable. Latex-free.

## **Order Information**

Lay hand flat and measure around the MP joints (large knuckles).

	Color	Right	Left	Fits (circumference at MP's)
Sm/Med	Black	P2006-R23BK	P2006-L23BK	6" to 7 1/4"
Sm/Med	Grey	P2006-R23GR	P2006-L23GR	6" to 7 1/4"
Med/Lrg	Black	P2006-R34BK	P2006-L34BK	7 1/4" to 9 1/4"
Med/Lrg	Grey	P2006-R34GR	P2006-L34GR	7 1/4" to 9 1/4"

Watch our **instructional video on the ThumSling** on our website:
www.3pointproducts.com/productvideos-of-splints/3pp-thumsling-video





Call 888-378-7763

Mon-Fri 8:30am - 5:00pm EST for ordering details

www.3pointproducts.com