

BEDCANE

BECAUSE YOU DESERVE SAFETY AND STYLE



PRODUCT #2041















FEATURES & BENEFITS

- Ergonomic No-Slip Handle Allows for easy transfer in and out of bed.
- **4-Pocket Organizer Included** Prevents entrapment and provides storage space for handy items.
- Portable Folds for travel or storage.
- Safety Strap Secures rail to any home or hospital bed.

SPECIFICATIONS

Rail Height: Adjustable 19" - 22"

Rail Width: 15"

Base Dimensions: 22" x 16" Weight of Product: 12 lbs Weight Capacity: 350 lbs

Made for Bed Type: Attaches to any home or hospital

bed with included safety strap

Assembly: Installs in minutes with 3 bolts and safety pin

Package Dimensions: $16.5" \times 24" \times 2.5"$; 10 lbs Case of 4 Dimensions: $25" \times 17.5" \times 11.5"$; 42.5 lbs

EASY ASSEMBLY (No Additional Tools Required)



Attach Handle to Bracket with Bolt and Safety Lock Pin



Attach Pocket Organizer and Safety Strap



Secure to Bed Frame with Included Safety Strap