

AMFAB™

 AbilityOne

1446 South 35th Street
Galesburg, MI 49053
800.432.6322



Automatic Overbed Table/Model 1005

Long-Term Care Tables



AmFab's first automatic, low overbed table is designed for use with low beds with a 3-1/4" floor to frame clearance. Available in both the "H" style and the popular "U" style bases.

The 1005 features a standard 15" x 30" walnut melamine top that can be raised with the slightest upward pressure, and lowers with a touch of the convenient trip lever. AmFab's patented automatic spring-assisted lift mechanism allows infinite height adjustment between the 19-3/4" to 28" height range.

The frame and column are made of 16 and 17 gauge welded tubular steel with an attractive chrome finish for years of durable service. 40 mm dual wheel casters are appropriate for use on any surface for easy maneuverability.

The 1005U features a 24" open area between the legs of the base with the casters mounted at the extreme points for maximum stability.

The "U" base model also makes the perfect companion for activities including Hemi, Super-Hemi, and Pediatric height wheelchairs.

Height
19-3/4" to 28"

Weight
26 lbs., Cube 1.1, UPS Shippable

Shipping
One table per carton, knocked down. Table easily assembled with four bolts.

Options (Available at additional cost)

- Alternate high pressure laminate
- Vanity and mirror (factory installed when ordered with table)
- Painted frame*

* Please allow six weeks additional for delivery.

AMFAB™

a AbilityOne

1446 South 35th Street
Galesburg, MI 49053
800.432.6322