

AIR MATTRESSES



WE CARE LIKE FAMILY™

Reference Chart



Bariatric HD Airfloat



Airfloat 100



Airfloat 300



Airfloat 500



Airfloat Deluxe



Bariatric HD Airfloat Plus

Unique Features	Mattress Type					
	Bariatric HD Airfloat 10439	Airfloat 100 10441, 10448, 10449	Airfloat 300 10442	Airfloat 500 10447	Airfloat Deluxe 10443-10444	Bariatric HD Airfloat Plus 10445-10446
36 x 80	—	Yes	Yes	Yes	Yes	—
36 x 84 (Optional)	—	Yes	—	—	—	—
42 x 80	—	—	—	—	—	Yes
48 x 80	—	—	—	—	—	Yes
54 x 80	Yes	—	—	—	—	—
Raised Perimeter (Optional)	—	—	—	—	Yes	—
Single Layer Cell w/ Foam Base (Optional)	—	Yes	—	—	—	—
Cells	Cell-On-Cell	Single Cell	Single Cell	Single Cell	Cell-On-Cell	Cell-On-Cell
Cell Height	10"	8"	8"	8"	8"	10"
Pillow Function	Yes	Yes	Yes	Yes	Yes	Yes
CPR Valve	Yes	Yes	Yes	Yes	Yes	Yes
Cover type	Nylon Quilted	Nylon	Nylon Quilted	Nylon	Nylon Quilted	Nylon Quilted
Manual Dial Pump	—	Yes	Yes	—	—	—
Digital Pump	Yes	—	—	Yes	Yes	Yes
Alternating Mode	Yes	Yes	Yes	Yes	Yes	Yes
Static Mode	Yes	Yes	Yes	Yes	Yes	Yes
Alarm Type	Audible and Visual	Visual	Audible	Audible	Audible	Audible
Max Firm Function	Yes	—	—	—	Yes	Yes
Cycle Time (Minutes)	10, 15, 20, 25	10 Min. Fixed	10 Min. Fixed	10 Min. Fixed	10, 15, 20, 25	10, 15, 20, 25
Lock Button for Weight Setting	Yes	—	—	Yes	Yes	Yes
Seat Inflate (Fowler position)	Yes	—	—	Yes	Yes	Yes
Warranty	36 Months	1 Year Limited	1 Year Limited	18 Months Limited	24 Months Limited	24 Months Limited
Weight Capacity	1000 lb.	400 lb.	400 lb.	350 lb.	400 lb.	880 lb.