

# Latex Powder Free PL Examination Glove

Our OmniTrust Latex Powder Free PL is a reduced protein glove. It has our famous (micro-texture) matte finish for added grip control without loss of tactile sensitivity. This glove is manufactured with a non-chlorination process thus maintaining the stretch characteristics of latex. A unique polymer coating gives this glove excellent donning characteristics. Omni International latex gloves are flat packed for optimum dispensing in sizes small through extra large.



- Non-Chlorinated Latex
- Polymer Coated for Easy Donning
- Ambidextrous, Beaded Cuff
- Micro-Textured for Superior Grip Control
- Single Use Only

| Tensile Strength - Aged/Unaged | Elongation - Aged/Unaged |
|--------------------------------|--------------------------|
| min. 14 MPa / min. 18 MPa      | min. 500% / min. 650%    |

| Dimensions         | Finger                  | Palm                    | Cuff                    | Length      | Width      |
|--------------------|-------------------------|-------------------------|-------------------------|-------------|------------|
| <b>Small</b>       | min 0.08 mm<br>3.15 mil | min 0.08 mm<br>3.15 mil | min 0.08 mm<br>3.15 mil | min. 240 mm | 80 ± 5 mm  |
| <b>Medium</b>      | min 0.08 mm<br>3.15 mil | min 0.08 mm<br>3.15 mil | min 0.08 mm<br>3.15 mil | min. 240 mm | 95 ± 5 mm  |
| <b>Large</b>       | min 0.08 mm<br>3.15 mil | min 0.08 mm<br>3.15 mil | min 0.08 mm<br>3.15 mil | min. 240 mm | 106 ± 5 mm |
| <b>Extra Large</b> | min 0.08 mm<br>3.15 mil | min 0.08 mm<br>3.15 mil | min 0.08 mm<br>3.15 mil | min. 240 mm | 116 ± 5 mm |

| Item Numbers      | Small  | Medium | Large  | X-Large |
|-------------------|--------|--------|--------|---------|
| <b>113 Series</b> | 113-01 | 113-02 | 113-03 | 113-04  |

| Packaging          | Boxes      | Cases    |
|--------------------|------------|----------|
| packaged by weight | 100 Gloves | 10 Boxes |

