

Nitrile Exam Gloves Sterile Pairs Extended Cuff Series 114

Nitrile synthetic contains no natural rubber latex proteins. Soft, highly elastic formulation provides the maximum in comfort. Pair packaged, inner wrapper in ozone resistant multi-vac packaging.

Characteristics: Sterile, Powder-Free, Blue

SIZE	ITEM NO.
XS	114050
S	114100
М	114200
L	114300
XL	114350

All Sizes: 50 Pairs/Box, 4 Boxes/Case

Tested for use with Chemotherapy Drugs



Features:

- Textured Finger = an Improved Wet/Dry Grip
- Low Modulus = a Softer, More Confortable Fit

Average Palm Thickness

Average Finger Thickness

- Non-Latex = No Risk of Latex Allegens
- Extended Cuff = Better Coverage

Manufactured in compliance with multiple international standards, including ASTM D6319, ASTM D5151, and ASTM F1671

Tested for use with Chemotherapy Drugs as per ASTM D6978 Standard Practice for Assessment of Medical Gloves to Permeation by Chemotherapy Drugs <i>Minimum Breakthrough Detection Time in Minutes</i> :		Tensile Strength & Elongation	Before Aging	After Accelerated Aging		
Ŭ		E 11 (50.0 (1)	. 0.10	Tensile Strength (Mpa)	32	33
Carmustine (BCNU) (3.3 mg/ml) Cisplatin (1.0 mg/ml)	43.8 >240	Fluorouracil (50.0 mg/ml) Methotrexate (25.0 mg/ml)	>240 >240	ASTM Requirement Min. (Mpa)	14	14
Cyclophosphamide (Cytoxan) (20.0 mg/ml)	>240	Mitomycin C (0.5 mg/ml) Paclitaxel (Taxol) (6.0 mg/ml)	>240 >240 97.8	Elongation (%)	691	667
Dacarbazine (10.0 mg/ml) Doxorubicin Hydrochloride (2.0 mg/ml)	>240 >240			ASTM Requirement Min. (%)	500	400
Etoposide (Toposar) (20.0 mg/ml)	>240	Vincristine Sulfate (1.0 mg/ml)	>240			
CAUTION: Please note Carmustine and Thio		ave extremely low permeation times	of 43.8	Average Length 11.5 in +	290 mm	

minutes and 97.8 minutes.

Caution: In compliance with ASTM D6978, the testing conditions used are intended to approximate the worst-case conditions for clinical use. Testing was conducted on single layer glove material. It is the users' responsibility to determine the applicability of these gloves for their intended use with chemotherapy drugs.

50 13485:2075 CERTIFICATION



Innovative Healthcare Corporation is certified to ISO 13485:2016 QMS for medical devices.



View this product on our website:

4.7 mil 🔶 0.12 mm

6.3 mil + 0.16 mm

Product Solutions You Trust