

## **Electric Youth Bed**

Model # 115E-H-A



## Features:

- Electric bed utilizes hand control for positioning
- Full safety length bedside rails that meet 2%" spacing requirments for preventing entrapment
- 5" braking casters
- Solid pan sleeping surface (non-helical design)
- Includes 4 IV hook holes and 4 urine bag hooks
- 5" Perforated Comfort Mattress (NOT INCLUDED)

## **Specifications**

■ Overall Dimensions: 72.5"L x 36"W

■ Height Adjustment: 23" - 33"

■ Sleeping Surface: 33.5" x 64.5"

■ Weight Capacity: 350 lbs.

Phone | 716.759.7200 Toll Free | 800.274.2742 Fax | 716.759.0700 Email | info@novummed.com www. novummed.com